**Instructions for Use**  
**RFA-GP-BAY**  
Adult Dispersive Electrode

Baylis Medical Company Inc.  
5959 Trans-Canada Highway  
Montreal, Quebec, Canada, H4T 1A1  
Tel: (514) 488-9801/ (800) 850-9801 Fax: (514) 488-7209  
www.baylismedical.com

**EU Authorized Representative:**  
Quality First International  
20 Eversley Road, Bexhill-on-Sea, East Sussex  
TN40 1HE, United Kingdom  
Tel: +44-(20)-8-522-1937 Fax: +44-(20)-8-522-1937

© Copyright Baylis Medical Company Inc, 2012-2014

The Baylis Medical logo is a trademark and/or a registered trademark of Baylis Medical Company Inc. in the United States of America and/or other countries.

**English**

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

**I. DEVICE DESCRIPTION**

The RFA-GP-BAY Adult Dispersive Electrode (Dispersive Electrode) is used with BMC Radiofrequency Puncture Generator (Generator) during RF procedures to provide a return connection from the patient to the generator

**II. WARNINGS**

- No modification of electrode is allowed. Modification may result in patient or operator harm.
- Failure to achieve good skin contact by the entire adhesive surface may result in an electrosurgical burn or poor electrosurgical performance.
- Do not coil dispersive electrode cable or allow the cable to overlie other electrosurgical monitoring cables or equipment as the unintended transfer of potentially harmful RF energy may result.
- Conductive parts of the electrode and associated connectors should not contact any other conductive parts including earth, as this contact may increase the risk of burns or harm to the patient or operator.
- Difficulty in achieving cutting energies requires evaluation. **Stop.** Do not proceed or increase power settings until you have checked all components of the electrosurgical circuit including the active electrode and its cable, the patient, dispersive electrode adherence and integrity and the Generator and its connectors. Indiscriminate power increase may result in patient burns.

**III. CAUTION**

- Heat applied by thermal blankets or other sources are cumulative with heat produced at the dispersive electrode. Choice of application site removed from other heat sources reduces the risk of patient injury.
- Always use the lowest power settings to achieved desired surgical effect.
- Do not use a needle monitoring electrode with dispersive electrodes.
- Avoid skin to skin contact when positioning patient, using dry gauze where necessary.
- When wrapping electrode around a limb, the electrode must not touch or overlap itself.
- Additional electrode gel is not required and should not be used.
- If patient is repositioned, re-inspect dispersive electrode and all connections.
- **DO NOT** reuse or relocate the dispersive electrode after initial application.
- Immediately wipe up any spilled fluids in area of the dispersive electrode. Keep application site dry.
- Avoid applying pressure to dispersive electrode through use of straps, tie downs, tape, etc.

**IV. DIRECTIONS FOR USE**

1. Select a well vascularized site in close proximity to the surgical site (lower back or thigh is recommended). Avoid placement on bony prominence, skin lesions, folds, scars, metal prosthesis or near EKG electrodes and cables. Do not apply where fluid may pool (e.g. under the buttocks).
2. Shave application site as required, clean to remove oil, lotions etc., and allow to dry thoroughly.
3. Do not open package until ready to apply electrode to skin. Inspect electrode and cable. Do not use if you suspect product is expired or unsafe (check expiration date printed on packaging).
4. After the patient is in the final position, carefully separate the electrode from the liner starting at the electrode edge farthest from the cable or cable attachment.

Avoid excess skin or finger contact with the adhesive surface prior to application. Electrode gel should not be used.

5. Apply electrode firmly, ensuring full adhesion and contact with the skin.
6. Connect the cable to the return electrode site on the generator.
7. To remove, gently support the skin and slowly peel away the electrode. Rapid removal may cause skin irritations.
8. Electrode should only be applied so that it is perpendicular to the direction of current flow from the operative site.
9. After use, this product may be a potential biohazard. Handle and dispose of in accordance with hospital or facility medical practice and following local, state and federal laws.

**V. LABELING AND SYMBOLS**

	Manufacturer		Do not reuse
	Latex Free		Lot Number
	Quantity per box		EU Authorized Representative
	Use By		Consult Instructions for Use
	Pouch- Tear open		Catalogue (Model) Number
	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.		
	Operating and storage temperature between 5°C and 35°C Operating and storage humidity between 30%-75%		
	Transportation Temperature between -29°C and 50°C Transportation humidity between 35%-85%		
	Shaded figure indicates suggested weight guidelines for Grounding Pads : >80lbs (>36.3 kg)		

**VI. LIMITED WARRANTY – Disposables and Accessories**

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date.

This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

**DISCLAIMER AND LIMITATION OF LIABILITY**

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE. THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product. This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products	The shelf life of the product
Accessory Products	90 days from the shipment date

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