BAYLIS MEDICAL COMPANY INC.

**BMC RADIOFREQUENCY PUNCTURE
GENERATOR**

EN Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

© Copyright Baylis Medical Company Inc., 2001-2012

PowerWire™, NRG™ and The Baylis Medical Company logo are trademarks and/or registered trademarks of Baylis Medical Company Inc. in the United States of America and/or other countries.

BAYLIS MEDICAL COMPANY INC.

BMC RADIOFREQUENCY PUNCTURE GENERATOR

Instructions for Use

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

Manufacturer:

Baylis Medical Company Inc
5959 Trans-Canada Highway
Montreal, Quebec
Canada
H4T 1A1

Phone: +1 (514) 488-9801
Fax: +1 (514) 488-7209

**EU Authorized
Representative:**

Quality First International
20 Eversley Road
Bexhill-on-Sea
East Sussex
TN40 1HE
United Kingdom

Phone: +44-(20)-8-522-1937
Fax: +44-(20)-8-522-1937

TABLE OF CONTENTS

TABLE OF CONTENTS	3
LIST OF FIGURES	5
LIST OF TABLES	5
PREFACE	6
SECTION 1: DEVICE DESCRIPTION	7
SECTION 2: INDICATIONS/CONTRAINDICATIONS	8
2.1. INDICATIONS FOR USE	8
2.2. CONTRAINDICATIONS	8
SECTION 3: WARNINGS, PRECAUTIONS, AND ADVERSE EVENTS	9
3.1. WARNINGS	9
3.2. PRECAUTIONS	10
3.3. ADVERSE EFFECTS	11
SECTION 4: UNPACKAGING AND REPACKAGING	12
4.1. UNPACKAGING	12
4.2. REPACKAGING	12
SECTION 5: CONTROLS AND DISPLAYS	13
5.1. GENERATOR SOFTWARE FLOW CHART	15
5.2. RF ON/OFF BUTTON/LIGHT	16
5.3. POWER	17
5.4. TIME	18
5.5. IMPEDANCE DISPLAY	19
5.6. SET INDICATOR	19
5.7. MEASURED INDICATOR	19
5.8. FAULT INDICATOR	20
5.9. POWER SWITCH	20
5.10. VOLTAGE SELECTOR SWITCH	20
5.11. FOOTSWITCH	20
SECTION 6: PORTS/CONNECTORS	22
6.1. ISOLATED PATIENT CONNECTIONS	22
6.2. FOOTSWITCH CONNECTOR	22
6.3. POWER CONNECTOR	22
6.4. SERIAL PORT	22
6.5. EQUIPOTENTIAL GROUND CONNECTION:	23
SECTION 7: DIRECTIONS FOR USE	24
7.1. READ INSTRUCTIONS FOR USE	24
7.2. VERIFY VOLTAGE SELECTOR SETTING	24
7.3. CONNECT GENERATOR POWER CORD	24
7.4. CONNECT BMC RFP CONNECTOR CABLE	24
7.5. CONNECT THE DIP ELECTRODE	25
7.6. TURN THE GENERATOR “ON”	25

7.7.	SELECT POWER LEVEL.....	26
7.8.	SET TIME DURATION	26
7.9.	CONFIRM GENERATOR AND ACCESSORY SETTINGS AND STATE	26
7.10.	RF POWER DELIVERY	26
7.11.	TURNING RF POWER “OFF”	27
7.12.	RE-APPLY RF POWER	27
SECTION 8: SERVICE AND MAINTENANCE		28
8.1.	CLEANING	28
8.2.	FUSE REPLACEMENT	28
SECTION 9: SPECIFICATIONS.....		29
9.1.	GENERAL SPECIFICATIONS.....	29
9.2.	POWER DELIVERY	31
9.3.	OUTPUT POWER GRAPHS.....	32
9.4.	ERROR CODES.....	34
9.5.	IEC ELECTRICAL SAFETY AND EMC SPECIFICATIONS	35
9.6.	LABELING AND SYMBOLS.....	39
SECTION 10: LIMITED WARRANTIES AND DISCLAIMER.....		41

LIST OF FIGURES

	Page
Figure 5-1 Generator Front Panel	13
Figure 5-2 Generator Back Panel.....	14
Figure 5-3 - Software Flow Chart.....	15
Figure 9-1 Power vs Load at 25W Setting.....	32
Figure 9-2 Power vs Load at 12W Setting.....	32
Figure 9-3 Peak Voltage vs Power Setting	33
Figure 9-4 Power Output vs Power Setting	33

LIST OF TABLES

	Page
Table 9.5-1 IEC Electrical Safety Specifications	35
Table 9.5-2 IEC EMC Specifications (Emissions)	35
Table 9.5-3 IEC EMC Specifications (Immunity).....	36
Table 9.5-4 IEC Recommended Separation of RF Communication Equipment	38

PREFACE

For the user's convenience, the BMC Radiofrequency Puncture Generator will be referred to in this Operator's Manual as the "Generator". The Generator may be used with active electrodes that have been separately cleared for use with the Generator. These separately cleared active electrode devices include, but are not limited to, the Nykanen Radiofrequency Wire, the *PowerWire™* Radiofrequency Guidewire, the Toronto Transseptal Catheter, and the *NRG™* Transseptal Needle - they will generally be referred to in this Operator's Manual as the "Catheter". The catheter is connected to the Generator through the appropriate BMC connector cable. The footswitch is an accessory to the BMC Radiofrequency Puncture Generator.

The use of the Baylis Medical Radiofrequency Puncture Generator is fully described in this manual, including a description of the Generator, its controls, displays, and a sequence for its operation. In addition, other information of importance to the user is supplied. For specific instructions pertaining to the use of any one of the separately cleared catheters, please refer to the instructions for use for the respective catheters.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

SECTION 1: DEVICE DESCRIPTION

The Baylis Medical Company Radiofrequency Puncture Generator is a component of the Baylis Medical Company Radiofrequency Puncture System. The Generator is operated in conjunction with a BMC Catheter and a commercially available external Disposable Indifferent (Dispersive) Patch (DIP) Electrode that meets or exceeds ANSI/AAMI Standard HF-18 and/or ISO 60601-2-2. The BMC Catheter delivers power in a monopolar mode between its distal tip electrode and the indifferent patch electrode. Detailed information regarding the Catheter is contained in a separate manual that accompanies each Catheter.

The Generator produces continuous radiofrequency (RF) power output at a fixed frequency in the range of 455 kHz to 465 kHz in a monopolar mode. Connections for the Catheter and a patient return electrode are provided. Controls on the front panel allow the power level and the duration of the RF output to be set. In addition, controls on the front panel provide on/off control of the output. Measured values for output power, impedance, and elapsed time are displayed during RF output whilst a constant audible tone is also produced. The Generator has several built-in safety features, such as warning messages, an automatic shutoff (if the measured impedance falls below 100 ohms or exceeds 6000 ohms), and maximum voltage and current limits.

SECTION 2: INDICATIONS/CONTRAINDICATIONS

2.1. INDICATIONS FOR USE

The Baylis Medical Radiofrequency Puncture System is intended to create an atrial septal defect in the heart.

2.2. CONTRAINDICATIONS

Usage of the Radiofrequency System is not recommended for usage other than the indicated use.

SECTION 3: WARNINGS, PRECAUTIONS, AND ADVERSE EVENTS

3.1. WARNINGS

- DO NOT attempt to operate the Generator before thoroughly reading this User's Manual. It is vital that the operating instructions for the equipment be read, understood, and followed properly. For future reference, retain this User's Manual in a convenient, readily accessible place.
- The Generator is intended for use with approved Catheters and accessories only.
- Do not remove the cover of the Generator. Removal of the cover may result in injury and/or damage to the Generator.
- When the Generator is activated, conducted and radiated electrical fields may interfere with other medical equipment. Care should be taken to limit the effects that electromagnetic interference (EMI) produced by the Generator has on other equipment.
- Laboratory staff and patients can undergo significant x-ray exposure during RF Puncture procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure.
- Do not attempt to puncture with an initial power setting greater than that recommended by the Catheter Instructions For Use. The output power selected should be as low as possible for the intended purpose.
- Failure of the Generator could result in an unintended increase of output power.
- Place monitoring electrodes as far away from the surgical site as possible, to avoid burns or interference with other equipment. The use of needle monitoring electrodes (or other small area electrodes) during RF output is not recommended. In all cases, monitoring systems incorporating high frequency current limiting devices are recommended.
- Skin-to-skin contact (for example between the arms and body of the patient) should be avoided, for example by insertion of dry gauze).
- The Generator is not equipped with a Contact Quality Monitor, thus loss of safe contact between the dispersive and the patient will not result in an auditory alarm.
- During RF output, implated devices such as pacemakers may be affected. Qualified advice should be obtained as necessary, to minimize the risk from injury to due implanted device malfunction.

3.2. PRECAUTIONS

- The Generator is intended for use with approved Catheters, Wires and accessories only. Ensure that the rated accessory voltage is equal to or greater than the Generator's maximum output voltage.
- Ensure that the catheter connector cables and dispersive electrode cables are positioned in such a way that contact with the patient or other leads is avoided.
- It is recommended not to exceed the specified number of power applications per Catheter, as indicated within the Catheter's specific instructions for use.
- Do not attempt to use the Generator before thoroughly reading the Generator User's Manual.
- Only physicians thoroughly trained in RF Puncture techniques, in a fully equipped catheterization laboratory, should perform RF Puncture procedures.
- Adequate filtering must be used to allow continuous monitoring of the surface electrocardiogram (ECG) during RF power applications.
- Read and follow the manufacturer's instructions for use of the Disposable Indifferent (Dispersive) Patch (DIP) electrode. Always use dispersive electrodes that meet or exceed IEC 60601-2-2 requirements. The entire area of the dispersive electrode should be reliably attached to the patient's body and as close to the operating field as possible.
- Placement of the dispersive electrode on the thigh could be associated with higher impedance.
- The Generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the Catheter and dispersive electrode, particularly when operating the device.
- During power delivery, the patient should not be allowed to come in contact with grounded metal surfaces or metal surfaces which have an appreciable capacitance to earth (for example operating table supports, etc). The use of antistatic sheeting is recommended for this purpose.
- Apparent low power output or failure of the equipment to function properly at normal settings may indicate faulty application of the dispersive electrode or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication.

- In order to prevent the risk of ignition make sure that flammable material is not present in the room during RF power application.
- Take precautions to limit the effects that the electromagnetic interference (EMI) produced by the Generator may have on the performance of other equipment.
- Regularly inspect and test re-usable cables and accessories.
- Baylis Medical Company relies on the physician to determine, assess, and communicate to each individual patient all foreseeable risks of the Generator.
- The mains power cord of the Generator must be connected to a properly grounded receptacle. Extension cords and/or adapter plugs must not be used.

3.3. ADVERSE EFFECTS

Adverse events that may occur while using the Baylis Medical Radiofrequency Puncture System include:

- Atrial Fibrillation and/or Atrial Flutter
- Myocardial Infarction
- Sustained arrhythmias leading to Ventricular Tachycardia.

SECTION 4: UNPACKAGING AND REPACKAGING

4.1. UNPACKAGING

The Generator's shipping carton contains all of the components identified below. Unpack the Generator and accessories carefully and visually inspect for damage. Notify the carrier immediately if the shipment carton is damaged. Verify that the following items are received:

- 1 Generator
- 1 Operator's Manual
- 1 Power Cord

Read the Directions for use in Section 7 of this manual very carefully and thoroughly. If there are any discrepancies or concerns, notify Baylis Medical Company. Store the shipping carton in a safe place for future use.

4.2. REPACKAGING

If it is necessary to repack and ship the Generator, use the original shipping carton and packing materials to ensure that no breakage occurs. Disconnect all the cables and accessories and place them into the locations in the carton that are reserved for these units. Do not attempt to forcefully fit all the components in the carton.

SECTION 5: CONTROLS AND DISPLAYS

The generator has controls and displays on both its front and rear panels as shown in Figures 5-1 and 5-2, respectively. Related controls and displays are grouped into distinct areas and are discussed below.


Figure 5-1 Generator Front Panel


Figure 5-2 Generator Back Panel

5.1. GENERATOR SOFTWARE FLOW CHART


Figure 5-3 - Software Flow Chart

FRONT PANEL

5.2. RF ON/OFF BUTTON/LIGHT

The RF ON/OFF Button/Light is shown in Figure 5-1G. This button is responsible for the three operational modes of the Generator: “READY”, “ON”, and “DONE”. A relationship of the modes and the events required to select each mode is provided in section 5.1

In the “READY” mode the power level and the count-up timer can be set, the impedance display is blank, the RF ON/OFF Button is not illuminated, and the audible tone is off.

In the “ON” mode, which is initiated when the RF ON/OFF Button is pressed, RF output is active, RF ON/OFF Button/Light is illuminated green, timer display counts up from 0 seconds to the set time, the audible tone is ON, and the actual power and impedance are displayed. The RF output is terminated and the mode changes to DONE (or FAULT) mode when the timer elapses, when the RF ON/OFF Button is depressed during ON mode, or when an error is detected.

In the DONE mode, the RF output is OFF, the RF ON/OFF Button is not illuminated, and the audible tone is off. The Generator passes through the DONE mode as a transition from the ON mode back to the READY mode. The DONE mode lasts for five seconds and then automatically changes to the READY mode. The last measured Power, Impedance, and Timer Displays are held during the DONE mode. The last two seconds of the DONE mode will display an error code if termination was due to an error.

Note: Two other modes that the Generator has are the POST mode and the FAULT mode.

The Generator enters into the POST mode when it is first turned on. In the POST mode the Generator performs a self-test of power generation, measurement, and control circuitry. The Generator either transitions to the READY mode upon successful completion of a self-test or transitions to the FAULT mode if the self-test fails.

The Generator enters into the FAULT mode when a non-recoverable error is detected, at which point the FAULT light (Fig. 5-1H) is illuminated and a rapid burst of audible tones is produced.

5.3. POWER

- **POWER DISPLAY**

The POWER Display is shown in Figure 5-1A. The power output is displayed in 1-Watt increments and the range is 0-25 Watts. The POWER Display is a two digit, green LED numeric display. When the generator is initially powered “ON”, the green LEDs will display the default power set point of 0 Watts. The desired power output can be set using the POWER ▲/▼ Buttons (Fig. 5-1B/5-1C). The POWER Display shows the set power (1-25 Watts) during the READY mode and then displays the measured delivered power (0-25 Watts) during ON mode. The measured power may be lower than the set power, depending on tissue impedance. The POWER Display will show the last measured power (0-25 Watts) during the DONE mode for five seconds while the generator completes the DONE mode and transitions to the READY mode. When the generator enters READY mode, the prior preset power level will be displayed, which can be adjusted as desired.

- **POWER ▲/▼ BUTTONS**

The POWER ▲/▼ Buttons are shown in Fig. 5-1B/5-1C. By pressing the POWER ▲/▼ Buttons you can set the level of power you want the generator to deliver to the Catheter during the Puncture process. By pressing the button once, you increase/decrease the power setting by 1 Watt. If you continue to hold down the button, you will see the power level scroll up/down rapidly. Once the button is released, the scrolling will stop. The POWER Display (Fig. 5-1A) shows the maximum amount of power the generator will deliver to the Catheter. The range for the POWER is 1-25 Watts in increments of 1 Watt.

Note: During the process of application of energy, the POWER ▲/▼ Buttons can be pressed to increase/decrease the level of power delivered by 1-Watt increments. The POWER Display will show the actual power delivered during the application of radiofrequency energy. The scrolling function of the POWER ▲/ ▼ Button will not work during the delivery of RF power (RF Power “ON” mode).

5.4. TIME

- **TIME DISPLAY**

The TIME Display is shown in Figure 5-1E. This display is a two-digit, green LED numeric display. The TIME Display will read zero when the Generator is initially powered “ON”. The desired output duration expressed in seconds can be set by the TIME ▲/▼ Buttons (Fig. 5-1D/5-1F), during the READY mode. The count-up time duration (0-99 seconds) is shown during the ON mode, as RF energy is delivered. If RF power delivery is stopped before the preset duration has been reached, the actual duration of RF power delivery will be displayed for five seconds while the generator reaches the DONE mode and transitions to the READY mode. Once in READY mode, the prior preset time duration will be displayed and can be adjusted as desired.

- **TIME ▲/▼ BUTTONS**

The TIME ▲/▼ Buttons are shown in Figure 5-1D/5-1F. Pressing the ▲/▼ buttons increase/decrease the duration that RF energy will be delivered. By pressing the button once, you can increase the maximum duration by incremental steps of 1 second. By continuously holding down the button, you will see the numerical display scroll rapidly.

Note: The TIME ▲/ ▼ Buttons are NOT functional during the RF Power “ON” mode (when RF power is being delivered).

5.5. IMPEDANCE DISPLAY

The IMPEDANCE Display is shown in Figure 5-1K. The IMPEDANCE Display is a four-digit, green LED numeric display. The impedance measurement range is 100-6000 ohms. The displayed value will fluctuate during impedance measurement since it is a measure of catheter-to-tissue contact. During the READY mode, the IMPEDANCE Display remains blank since impedance is only measured during the ON mode while RF output is active. When RF power is discontinued (DONE mode) the last measured value of the impedance will remain on the display for five seconds while the generator completes the DONE mode and transitions to the READY mode. Once in the READY mode, the Impedance Display becomes blank.

Note: If the impedance exceeds 6000 ohms, the Generator will automatically cease the delivery of RF power and the IMPEDANCE Display will show a flashing “E002” error message to alert the user to the high impedance value. This error message will also be seen if the user attempts to deliver RF power to the Catheter when there is an electrical open circuit. The open circuit may be caused by various influences, such as a Catheter that was not connected to the BMC Connector Cable. If the impedance drops below 100 ohms, the Generator will automatically cease the delivery of RF power and the display will show a flashing “E001” error message to alert the user to the low impedance value.

5.6. SET INDICATOR

The SET Indicator is shown in Figure 5-1M. It is a yellow LED that is illuminated during the READY mode. When the generator is in the READY mode, the timer and power levels can be set.

5.7. MEASURED INDICATOR

The MEASURED Indicator is shown in Figure 5-1N. It is a green LED that is illuminated when the Generator is in the ON mode. The Measured Indicator switches off after the Generator transitions from the DONE mode to the READY mode, but stays lit during the DONE mode.

5.8. FAULT INDICATOR

The FAULT Indicator is shown in Figure 5-1H. It is a red LED that indicates that a system error has occurred. The error code is displayed on the IMPEDANCE Display (Fig. 5-1K). There are three categories of error codes: error(s) during the POST mode, error(s) during the ON mode, and error(s) in the CPU Hardware self-test. Recoverable errors do not cause the Fault Indicator to illuminate, but instead cause a transition to the DONE mode. Non-recoverable errors cause a transition to the FAULT mode and illuminate the LED. To interpret the error codes, refer to the table entitled “Summary of Error Conditions and Error Codes” in Section 9.4. RF power is disabled when the FAULT Indicator is illuminated.

5.9. POWER SWITCH

The POWER Switch is shown in Figure 5-2A. Turn the Generator “ON” by pressing the rocker switch to the “I” position. Pressing the rocker switch to the “0” position will power the Generator “OFF”.

5.10. VOLTAGE SELECTOR SWITCH

The orientation of this switch (Figure 5-2D), in conjunction with the Fuse Drawer (Figure 5-2E), determines the AC input voltage range (100-120V~ or 220-240V~).

Note: The voltage selector and the orientation of the Fuse Drawer must be set to the same voltage setting. Do not change the position of the voltage selector while the system is plugged in.

5.11. FOOTSWITCH

The pneumatic footswitch can be used as an alternative to the RF ON/OFF button described in section 5.2. As an alternative, it is also responsible for the three operational modes of the Generator: “READY”, “ON”, and “DONE”. A relationship of the modes and the events required to select each mode is provided in section 5.1. Footswitch operation is abbreviated on the label “Generator Operation” found on the top of the generator and described below.


In the “READY” mode the power level and the count-up timer can be set, the impedance display is blank, the RF ON/OFF Button is not illuminated, and the audible tone is off.

In the “ON” mode, which is initiated when the footswitch is pressed and held, RF output is active, RF ON/OFF Button/Light is illuminated green, timer display counts up from 0 seconds to the set time, the audible tone is ON, and the actual power and impedance are displayed. The RF output is terminated and the mode changes to DONE (or FAULT) mode when the timer elapses, when the footswitch

is released, when the RF ON/OFF Button is depressed during ON mode, or when an error is detected.

In the “DONE” mode, the RF output is OFF, the RF ON/OFF Button is not illuminated, and the audible tone is off. The Generator passes through the DONE mode as a transition from the ON mode back to the READY mode. The DONE mode lasts for five seconds and then automatically changes to the READY mode. The last measured Power, Impedance, and Timer Displays are held during the DONE mode. The last two seconds of the DONE mode will display an error code if termination was due to an error.

- **FOOTSWITCH LABEL**

SECTION 6: PORTS/CONNECTORS

6.1. ISOLATED PATIENT CONNECTIONS

CATHETER ISOLATED PATIENT CONNECTOR

The CATHETER ISOLATED PATIENT CONNECTOR is shown in Figure 5-1I. This port provides for the connection of the BMC Connector Cable to the Generator. The Connector Cable can be attached to the Catheter as specified in the Instructions for Use that is packaged with the Catheter. The Connector Cable MUST be installed in the CATHETER ISOLATED PATIENT CONNECTOR before the Generator will deliver any RF power. If the Generator is switched on without the Connector Cable connected, RF output will not start and the Impedance Display (Fig. 5-1K) will read E002. (See Error Codes – Section 9.4)

RETURN ELECTRODE ISOLATED PATIENT CONNECTOR

The RETURN ELECTRODE ISOLATED PATIENT CONNECTOR is shown in Figure 5-1J. This port provides for the connection of the Disposable Indifferent (Dispersive) Patch (DIP) electrode to the Generator. The DIP Electrode Cable uses a flat connector that is keyed for proper alignment. The DIP Electrode Cable MUST be installed in the RETURN ELECTRODE ISOLATED PATIENT CONNECTOR before the Generator will deliver any RF power. If the Generator is switched on without the DIP electrode Cable connected, RF output will not start and the Impedance Display (Fig. 5-1K) will read E002. (See Error Codes – Section 9.4)

6.2. FOOTSWITCH CONNECTOR

The FOOTSWITCH CONNECTOR is shown in Fig. 5-1L. This port provides for the connection of a footswitch to the Generator.

6.3. POWER CONNECTOR

The Power Connector is shown in Figure 5-2B it is a Hospital Grade connector.

6.4. SERIAL PORT

The Serial Port is shown in Figure 5-2C and is labeled “RS-232”. It is provided for computer interface to Baylis Medical Company approved accessory computers only.

6.5. EQUIPOTENTIAL GROUND CONNECTION:

This connector (Figure 5-2F) is attached to the chassis/earth ground. It is intended for earth reference connection in environments where equipotential ground cabling is used.

SECTION 7: DIRECTIONS FOR USE

7.1. READ INSTRUCTIONS FOR USE

Do not operate the generator or accessories before thoroughly reading their respective Instructions For Use. It is vital that the Instructions For Use for the equipment is read, understood, and followed properly.

7.2. VERIFY VOLTAGE SELECTOR SETTING

Verify that the voltage selector setting and fuse draw orientation match the supply voltage.

7.3. CONNECT GENERATOR POWER CORD

Connect the Generator Power Cord plug into a properly grounded AC electrical outlet. To ensure proper grounding, the Power Cord plug must be installed in an AC electrical wall outlet designated “Hospital Grade” or “Hospital Only”. **Never use an outlet without a grounding connection.** Position the generator for easy viewing of its front panel displays.

7.4. CONNECT BMC RFP CONNECTOR CABLE

Install the BMC RFP Connector Cable into the “CATHETER” socket located within the “ISOLATED PATIENT CONNECTIONS” area on the Generator front panel. Gently line up the connector pins with the socket and push in until the connector fits firmly into the socket. Secure the connector in place by rotating the locking ring in a clockwise direction.

To disconnect the BMC RFP Connector Cable from the Generator, rotate the locking ring counter clockwise to unlock, then grasp the connector and gently pull it straight out of the socket.

Never disconnect the BMC RFP Connector Cable by pulling on the cable cord.

Do not twist the BMC RFP Connector Cable while inserting or removing it from the Generator socket, or else the pin connections may be damaged.

To connect the Catheter to the BMC RFP Connector Cable, please refer to the individual Instructions For Use of each device.

Note: Position the Generator near the table where the procedure is to be performed. The generator may only be connected to recording systems providing patient electrical isolation in accordance with EN-60601.

7.5. CONNECT THE DIP ELECTRODE

Connect the DIP Electrode connector to the “RETURN ELECTRODE” connector found within the “ISOLATED PATIENT CONNECTIONS” area on the Generator front panel. Gently push the DIP electrode-fitting straight in, until firmly in place. To unplug, grasp the DIP electrode fitting and gently pull it out from the receptacle.

Before use, it is important to check the DIP Electrode for damage to the sealed foil package. Exposure to air, due to a damaged package, could cause the DIP Electrode to become dry and limit its grounding capability. Be sure the pad is moist and sticky to the touch before placing it on the patient. Do not attempt to relocate the patient grounding pad after its initial application. Electrode gel is NOT required and should NOT be used.

The DIP Electrode should be placed on a well-vascularized, convex skin surface. Do not place this electrode on the thigh, since this location is associated with higher impedance. Avoid scar tissue, bony prominence, adipose tissue, or any areas where fluid may pool. Shave, clean, and dry the application site as needed. Check for wrinkles or folds when applying the pad, as these will decrease conductivity.

7.6. TURN THE GENERATOR “ON”

Turn the Generator “ON” by pressing the rocker switch located on the Generator back panel (See Fig. 5-2A) to the “I” position. The Generator first performs a self-test of power generation, measurement, and control circuitry, as indicated by a constant audible tone. If all is satisfactory, the Generator enters into the READY mode, as indicated by the “SET” indicator being lit. If the self-test fails, the Generator enters into the FAULT mode, the red “FAULT” indicator will be illuminated, in addition, a rapid burst of audible tones will be emitted.

If the FAULT mode is entered (i.e. a system malfunction is detected during self-test), the Generator will not operate. To clear any malfunctions found during the self-test, the Generator must be powered “OFF” and then back “ON” with the self-test repeated. If the Generator fails again, the Generator cannot be used properly and Baylis Medical Company should be contacted for service. The Generator will NOT operate unless the initial self-test has been successfully completed.

Data Capture Option:

If a computer and recording software (an option available through Baylis Medical) are utilized, they will enable the user to capture the measured data generated during the RF puncture procedure. To enable the data capture option, connect the male end of the

RS-232 computer interface cable to the RS-232 port on the back of the generator (illustrated in Figure 5-2C) and connect the female end to the designated serial port on the computer. For instructions on loading the software and its operation sequence, please review the software's specific Instructions For Use.

7.7. SELECT POWER LEVEL

Set the desired RF power level (in Watts), as recommended by the Catheter Instructions for Use, by using the POWER ▲/▼ Buttons (Fig. 5-1B/5-1C), as described in section 5.3.

7.8. SET TIME DURATION

Set the desired or catheter Instructions For Use specified duration of the RF Power application (in seconds) by using the TIME ▲/▼ Buttons (Fig. 5-1D/5-1F). The TIME Display (Fig. 5-1E) will show the set duration. The duration of RF Power application can only be changed during the READY mode by using the TIME ▲/▼ Buttons, as described in section 5.4

Note: The TIME ▲/▼ Buttons are inoperative during RF power delivery.

7.9. CONFIRM GENERATOR AND ACCESSORY SETTINGS AND STATE

Prior to delivering RF power to the Catheter, first check that all the connections have been made properly and that the Power and Time settings are correct. Confirm that all of the requirements specified in the individual accessories' Instructions for Use have been met. Only once all the above conditions are met should one proceed to deliver RF power.

7.10. RF POWER DELIVERY

RF power is delivered when the RF ON/OFF Button (Fig. 5-1G) is pressed once. Pressing the RF ON/OFF Button again before the timer elapses will prematurely terminate RF power delivery. When RF power is delivered to the Catheter, the "POWER", "IMPEDANCE", AND "TIME" displays function as described in SECTION 5: and the RF ON/OFF BUTTON/LIGHT is illuminated green. It remains lit until the RF power delivery is discontinued by any of these occurrences:

- 1) Reaching the end of the set TIME duration.
- 2) Pressing the RF ON/OFF Button again.
- 3) When an error is detected.

If an excessively high impedance value is detected during delivery of RF power, a build-up of coagulum on the Catheter tip may be the cause and RF power output should be discontinued.

7.11. TURNING RF POWER “OFF”

RF output is terminated and the mode changes to the DONE mode when either the timer elapses or when the RF ON/OFF Button is pressed during the ON mode, at which point the RF ON/OFF Light will turn “OFF”.

RF output is also terminated and the mode changes to the FAULT mode when an error is detected. If the error is recoverable (E001 or E002), the unit proceeds to the DONE mode. For a non-recoverable error, which is illustrated by the illumination of the FAULT Indicator and a rapid burst of audible tones, the Generator must be turned “OFF” and then “ON” again to clear the error, after which the sequence of initializing steps must be performed again from Section 7.7

Note: If a non-recoverable error is encountered repeatedly, the Generator is not functioning properly and needs servicing. Contact Baylis Medical Company.

7.12. RE-APPLY RF POWER

To re-apply RF power, repeat steps from section 7.7

SECTION 8: SERVICE AND MAINTENANCE

The Generator requires no routine service or maintenance. If the Generator fails to operate when plugged into a proper AC power receptacle and the POWER Switch (Fig. 5-2A) is turned “ON”, a fuse may have blown. Contact Baylis Medical Company. The Generator contains no user-serviceable parts. Disassembly and attempted repair by unqualified personnel may create a hazardous condition and will void the warranty.

WARNING: DO NOT remove the cover of the Generator. Removing the cover may result in personnel injury and/or damage to the Generator.

8.1. CLEANING

The outer surface of the Generator may be cleaned with a mild soapy solution. DO NOT immerse the Generator or its accessories in any liquid. Avoid caustic or abrasive cleaners. If disinfecting is required, isopropyl alcohol may be used to clean the outer surfaces.

8.2. FUSE REPLACEMENT

1. Unplug the Power Cord from the Generator.
2. Use a precision slot screwdriver to remove the fuse drawer
3. Remove BOTH fuses from the fuse drawer and discard them.
4. Select TWO new fuses according to the specifications in Section 9.1, according to the supply voltage:
 - a. 100-120 VAC (Domestic Version – Canada/US):

<i>Current</i>	<i>Voltage Rated</i>	<i>Fuse Type</i>	<i>Melting I²t</i>	<i>Markings</i>
2.0A	250V	Slow Blow	14 I ² t	UL/CSA

- b. 220-240 VAC (International Version):

<i>Current</i>	<i>Voltage Rated</i>	<i>Fuse Type</i>	<i>Melting I²t</i>	<i>Markings</i>
1.0A	250V	Time lag	11 I ² t	IEC

5. Ensure the integrity of the new fuses by inspecting for physical damage that could affect the function of the fuse. Replace if either or both look damaged.
6. Replace the fuses in the fuse drawer.
7. Replace the fuse drawer back into the Generator in the orientation that **MATCHES THE SUPPLY VOLTAGE.**

WARNING: A mismatch of fuse type, Fuse Drawer orientation and Voltage Selector Switch can result in permanent damage to the Generator!

SECTION 9: SPECIFICATIONS

9.1. GENERAL SPECIFICATIONS

Line Power: 100-120 VAC (Domestic Version - CANADA/US)
220-240 VAC (International Version)

Current Rating: 1.0A, 60 Hz (100 - 120 VAC) - Domestic Version
0.5A, 50 Hz (220-240 VAC) - International Version

Fuse Rating: 2.0A, 250V, Slow Blow, 14 I²t, UL/ CSA (Domestic Version).
Recommended 2.0A fuse: Littlefuse, 218 series
(P/N: 0218002.HXP)
1.0A, 250V, Time Lag, 11 I²t, IEC (International Version)
Recommended 1.0A fuse: Littlefuse, 239 series
(P/N: 0239001.MXP)

Power Cord Length: 10 feet

Catheter Connector: Quick Connect female 4 pin

Indifferent Electrode Connector: Standard male 2-pin for commercial pads

Recommended Pads: Valley Lab, Pfizer

RS232 Cable Connector: Standard RS232 9-pin female Serial Port Connector

Generator:

Power Output: 1-25 Watts, settable in 1-Watt increments
25 Watts into a resistive load of 500-2000 ohms
Maximum output current of 0.224 A RMS
Maximum output voltage of 224.0 V RMS
Into resistive load range of 100-6000 ohms
Measurement accuracy: $\pm 5\% \pm 2$ Watts

Impedance: Measurement range 100-6000 ohms
Measurement accuracy: $\pm 10\% \pm 10$ ohms

Count-up Timer: Settable from 1-99 seconds
Display resolution: 1 second
Accuracy: 0.1 second

Dimensions: Width: 12 inches (30 cm)
 Length: 14 inches (35 cm)
 Height: 6 inches (15 cm)

Weight: 20 lb. (9.09 kg)

Leakage Current Measurements (No Fault Condition):

Catheter Source Current < 10 uA
 Indifferent Electrode Source Current < 10 uA
 Catheter Sink Current < 10 uA
 Indifferent Electrode Sink Current < 10 uA
 Case Leakage (Ground Open, Power Normal) < 100 uA
 Case Leakage (Ground Open, Power Reversed) < 100 uA

Dielectric Withstand (Hi-Pot) Test:

Mains to Catheter (4,000 V, 1min) PASS
 Mains to Indifferent Electrode (4,000 V, 1min) PASS

Accessory Rated Voltage:

Accessory Device	Rated Accessory Voltage (V_{rms})
<i>PowerWire</i> TM Radiofrequency Guidewire	224
Nykanen Radiofrequency Wire	224
<i>NRG</i> TM Transseptal Needle	224
Toronto Transseptal Catheter	224
Toronto Transseptal Wire	224

9.2. POWER DELIVERY

As an enhanced patient safety feature, both the output voltage and current of the Generator has been limited. Under normal circumstances (i.e. normal patient impedances), the maximum output power of the Generator is 25 Watts. In the event of unusually high impedance, the voltage limiting circuitry will limit maximum power. In the event of abnormally low tissue impedance, the current limiting circuitry will limit maximum power. Note that for most punctures, the voltage or current limits will not be reached. Also, note that ONLY the maximum power is affected. If the power setting is less than the maximum power available, the Generator will deliver the set power value.

Environmental Specifications

Storage

- Temperature: -20°C to 50°C. The unit should be gradually returned to the operating temperature range before use and stabilized for one hour before operation
- Relative Humidity: 10 % to 100 %, non-condensing
- Atmospheric Pressure: 500 to 1060 millibar

Operating

- Temperature: 15°C to 40°C
- Relative Humidity: 30 % to 75 %, non-condensing
- Atmospheric Pressure: 700 to 1060 millibar

9.3. OUTPUT POWER GRAPHS


Figure 9-1 Power vs Load at 25W Setting


Figure 9-2 Power vs Load at 12W Setting


Figure 9-3 Peak Voltage vs Power Setting


Figure 9-4 Power Output vs Power Setting

9.4. ERROR CODES

Summary of Error Conditions and Error Codes			
CPU Hardware self-test errors			
Error Code	Reason	Non-recoverable	Recoverable
H001 (I)	CPU configuration error	✓	
H002 (I)	COP time-out error	✓	
H003 (I)	ROM CRC error	✓	
H004 (I)	RAM error	✓	
H005 (I)	Illegal CPU instruction	✓	
POST errors			
P001 (I)	RF generation / measurement error	✓	
ON mode errors			
E001 (I)	Low impedance measured		✓
E002 (I)	High impedance measured		✓
E003 (P)	RF power too high	✓	
E004 (P)	RF current too high	✓	
E005 (P)	RF current measurement error	✓	
E006 (P)	RF voltage too high	✓	
E007 (P)	RF hardware shutdown	✓	

9.5. IEC ELECTRICAL SAFETY AND EMC SPECIFICATIONS

Table 9.5-1 IEC Electrical Safety Specifications

Device Description	
Class I, Defibrillation proof Type CF Equipment, IPX0, not AP/APG	
Mode of Operation: Continuous	
Electrical Isolation	<ul style="list-style-type: none"> Leakage current conforms to IEC 60601-1 Dielectric withstanding voltage conforms to IEC 60601-1
<p>EMC Emissions and Susceptibility: The BMC Radiofrequency Puncture Generator has been tested and found to comply with the limits for medical devices to the IEC 60601-1-2:2004. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This Generator generates, uses, and can radiate radiofrequency energy and, if not installed and used in accordance with the instruction given below, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation.</p>	

Table 9.5-2 IEC EMC Specifications (Emissions)

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The BMC Radiofrequency Puncture Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the BMC Radiofrequency Puncture Generator should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	The BMC Radiofrequency Puncture Generator must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class A	The BMC Radiofrequency Puncture Generator is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 9.5-3 IEC EMC Specifications (Immunity)

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The BMC Radiofrequency Puncture Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the BMC Radiofrequency Puncture Generator should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	<5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BMC Radiofrequency Puncture Generator requires continued operation during power mains interruptions, it is recommended that the BMC Radiofrequency Puncture Generator be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity (continued)

The BMC Radiofrequency Puncture Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the BMC Radiofrequency Puncture Generator should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the BMC Radiofrequency Puncture Generator including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> <p>$d = [1.17]\sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = [2.33]\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.


Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BMC Radiofrequency Puncture Generator or any of its components are used exceeds the applicable RF compliance level above, the BMC Radiofrequency Puncture Generator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating components or the entire BMC Radiofrequency Puncture Generator.


^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 9.5-4 IEC Recommended Separation of RF Communication Equipment

Recommended separation distances between portable and mobile RF communications equipment and the BMC Radiofrequency Puncture Generator			
The BMC Radiofrequency Puncture Generator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BMC Radiofrequency Puncture Generator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and BMC Radiofrequency Puncture Generator as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = [\frac{3.5}{V_1}] \sqrt{P}$	$d = [\frac{3.5}{E_1}] \sqrt{P}$	$d = [\frac{7}{E_1}] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
50	8.25	8.25	16.50
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

9.6. LABELING AND SYMBOLS

<i>FRONT PANEL</i>	
Defibrillator-proof, Patient Isolated connections	
High Frequency Isolated Patient Circuit	
Dispersive Return Electrode connection	
BMC Connector Cable connection	
Down	
Up	
Footswitch	
Output OFF	
Output ON	
Time	
IMPEDANCE (Ω)	IMPEDANCE
FAULT	FAULT
SET	SET
MEASURED	MEASURED

REAR PANEL	
Power OFF	
Power ON	
Alternating Current	
Attention: consult accompanying documents	
Dangerous Voltage	
Earth Ground	
Fuses	
VOLTAGE SELECTOR	VOLTAGE SELECTOR
Non-Ionizing Radiation	
Manufacturer	
Authorized Representative in the European Community	
Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.	
Consult Operating Instructions	
Catalogue Number	
Serial Number	
RS-232 Serial Connection	10101
Please contact your distributor or the device manufacturer for recycling of this device.	

SECTION 10: LIMITED WARRANTIES AND DISCLAIMER

LIMITED WARRANTY – RF Generators

Baylis Medical Company Inc. warrants the RF Generator and Foot Switch against defects in materials and workmanship to the registered owner at the time of purchase. All components of the RF Generator and Foot Switch are covered by the warranty as described below, except the connector cables, catheters, guide-wires, and accessories, which are covered in their own manuals and have their own warranties. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of the product. The length of the warranty is: (i) for the RF Generator, 1 year from shipment date, and (ii) for the Foot Switch, 90 days from shipment date.

This Limited Warranty applies only to new, original factory delivered products that have been used for their normal and intended uses. BMC's limited warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly installed, operated or maintained contrary to BMC's instructions. This warranty does not apply to any unit which has been subject to misuse, neglect, improper installation or that which has been altered, adjusted or tampered with by any person other than Baylis Medical authorized personnel.

If upon examination by authorized service personnel, it is determined that the malfunction is due to misuse or abuse, warranty provisions will not apply. An estimate of the cost of repair work will be given to the customer prior to servicing and repairing the unit.

The customer is responsible for returning the defective equipment to Baylis Medical at 5959 Trans-Canada, Montreal, Quebec, H4T 1A1 or to a specified address if different at his or her own expense. The customer shall obtain a return authorization number before shipping the unit back. Baylis Medical at its sole discretion can repair the unit or ship a new one. The units are to be shipped freight pre-paid for both the warranty period and out of warranty.

If, upon examination, it is determined that the fault had been caused by misuse or abnormal conditions of operation, the repairs will be billed to the customer as out-of-warranty repairs.

Instruments repaired under Baylis Medical standard repair program will be issued a thirty-day warranty against defects in both materials and workmanship, provided the original warranty period has passed. Instruments submitted due to defects in materials and workmanship during the thirty-day warranty period will be repaired at no charge to the customer.

DISCLAIMER AND LIMITATION OF LIABILITY

THE LIMITED WARRANTY SET FORTH HEREIN IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, REMEDIES, OBLIGATIONS AND LIABILITIES OF BAYLIS MEDICAL, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE OR PURPOSE. ANY WARRANTY OTHER THAN SET FORTH HEREIN IS EXPRESSLY DISCLAIMED.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. THESE PRODUCTS ARE BEING SOLD ONLY FOR THE PURPOSE DESCRIBED HEREIN, AND SUCH WARRANTY ONLY RUNS TO THE PURCHASER. IN NO EVENT SHALL BAYLIS MEDICAL BE LIABLE FOR ANY BREACH OF WARRANTY IN ANY AMOUNT EXCEEDING THE PURCHASE PRICE OF THE PRODUCT. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE

BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

The warranty periods for Baylis Medical products are as follows:

RF Generator Baylis #: RFP-100-115, RFP-100-230	1 year from the shipment date
Foot Switch Baylis #: PMA-FS	90 days from the shipment date

This page is intentionally left blank

This page is intentionally left blank

This page is intentionally left blank

This page is intentionally left blank

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Quality First International
20 Eversley Road
Bexhill-on-Sea
East Sussex
TN40 1HE
Royaume-Uni

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