



# VersaCross Radiofrequency System Reduces Time to Left Atrial Access versus Conventional Mechanical Needle

## HIGHLIGHTS

The study found LAAC sheath delivery with the **VersaCross™** RF Transseptal Solution was:

- ▶ Efficient: Transseptal puncture and LAAC sheath delivery on average in under 7 mins.
- ▶ Exchangeless: Faster LA access by combining a starter wire, RF transseptal device, and exchange rail in a 3-in-1 solution.
- ▶ Effortless: Controlled RF puncture with a single wire.

## INTRODUCTION

- ▶ Left atrial (LA) catheterization requires numerous device exchange steps, and has associated risks and safety concerns.
- ▶ The **VersaCross™** RF Transseptal Solution (Baylis Medical\*) enables vascular cannulation, transseptal puncture (TSP), and device exchange using a single RF-tipped pigtail wire.

## METHODS

- ▶ Consecutive series of left atrial appendage closure (LAAC) using WATCHMAN™ (Boston Scientific) or Amulet™ (Abbott) devices were retrospectively evaluated.
- ▶ Femoral access was obtained for inferoposterior TSP using two methods:

### Conventional approach (n=10):

- Requiring a starter wire, sharp mechanical needle (BRK-1™ Transseptal Needle, Abbott), fixed curve sheath (Swartz™ Transseptal Guiding Introducers, Abbott), and stiff exchange wire (Amplatz Super Stiff™, Boston Scientific or ProTrack™ Pigtail Wire, Baylis Medical)


### VersaCross™ RF Transseptal Solution (n=10):

- Comprised of the **VersaCross™** RF Wire, Sheath, and Dilator
- ▶ Efficiency was assessed in terms of time from femoral access to TSP, delivery of LAAC sheath in the LA, device release, overall procedure, and fluoroscopy use.
- ▶ Safety was assessed in terms of intraprocedural and in-hospital complications.

## RESULTS

- ▶ LAAC success was 100% using both methods, with no complications.
- ▶ Significant improvement in LA access times using **VersaCross™** RF Transseptal Solution vs. conventional method:
  - Shorter time to TSP [4.1±2.5 min vs. 8.4±4.0 min (p=0.009)]
  - Less time for LAAC delivery sheath into LA [6.7±2.4 min vs. 13.4±5.4 min (p=0.002; Figure 1)]
- ▶ Trend for overall procedural improvement using **VersaCross™** RF Transseptal Solution vs. conventional method:
  - Shorter time to device release [23.7±6.4 min vs. 31.2±10.0 min (p=0.062)]
  - Less fluoroscopy use [7.2±2.2 min vs. 11.4±5.9 min (p=0.061)]

Time for LAAC Sheath Delivery



**Figure 1.** LAAC sheath delivery is two times faster using the **VersaCross™** RF Transseptal Solution than the conventional workflow.†

## DISCUSSION & CONCLUSIONS

- ▶ **VersaCross™** RF Transseptal Solution combines a starter wire, transseptal needle, and exchange guidewire for faster LA access, and may improve overall procedural efficiency.

\* A wholly-owned subsidiary of Boston Scientific Corporation.

† From femoral access to LAAC sheath delivery using **VersaCross™** RF Transseptal Solution compared to a conventional mechanical needle and fixed curve sheath (Inohara et al., 2021).

All trademarks are property of their respective owners. Patents Pending and/or issued. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, Contraindications, Warnings, and Instructions For Use can be found in the product labelling supplied with each device or at [www.baylismedical.com](http://www.baylismedical.com).

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.



**Boston  
Scientific**  
Advancing science for life™

© 2023 Boston Scientific Corporation  
or its affiliates. All rights reserved.  
EP-1548701-AA