Instructions for Use

NorthStar™ Coronary Sinus Guidewire

ENGLISH 1

FRANÇAIS 5

 Manufactured for:
Baylis Medical Company Inc.
5959 Trans-Canada Highway
Montreal, Quebec, Canada, H4T 1A1
Tel: (514) 488-9801/ (800) 850-9801
Fax: (514) 488-7209
www.baylismedical.com

© Copyright Baylis Medical Company Inc., 2007-2020

The Baylis Medical logo and NorthStar are trademarks and/or registered trademarks of Baylis Medical Company Inc. in the United States of America and/or other countries. All other trademarks or registered trademarks are property of their respective owners.

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

CAUTION

This device should be used only by physicians trained in vascular intervention or in the placement of coronary venous leads.

Do not use if package is damaged. 

Do not re-sterilize  or re-use 

Do not wipe with alcohol.

See Instructions for Use 

Sterilized with ethylene oxide gas 

Use the guidewire prior to the "Use by date" on the package label, preceded by the  symbol

Refer to the instructions supplied with any diagnostic or interventional devices or with implantable coronary venous leads to be used in conjunction with the **NorthStar™** Coronary Sinus Guidewire for their intended uses, contraindications, and potential complications.

DESCRIPTION

The **NorthStar™** Coronary Sinus Guidewire is a steerable guidewire with a shapeable distal tip and as an option, a preshaped "J" tip is available. The **NorthStar™** Coronary Sinus Guidewire is available in three different shaft types, Extra Support (ES), Distal Support (DS) and Extra Distal Support (EDS).

CONTENTS: 1 Guidewire

STORAGE: Store in a cool, dark, dry place.

INDICATIONS FOR USE

The **NorthStar™** Coronary Sinus Guidewires are intended for use in the coronary and peripheral vasculature to facilitate the selective placement of diagnostic and interventional catheters

with compatible guidewire lumen and coronary venous leads.

CONTRAINDICATIONS

The **NorthStar™** Coronary Sinus Guidewires are not intended for use in the cerebral vasculature.
Patients judged not acceptable for vascular intervention.

WARNING

- A guidewire is a delicate instrument and must not be advanced, withdrawn, or torqued if resistance is met.
- Guidewire manipulations must always be observed under fluoroscopy. If the guidewire is removed and is to be re-inserted, it must be inspected for signs of damage (weakened or kinked segments) prior to re-introduction.
- Do not re-introduce if guidewire is weakened or kinked.
- Warning: For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

PRECAUTIONS

- Failure to follow the instructions may compromise guidewire performance and result in complications.
- Prior to use, confirm compatibility of guidewire outer diameter with the interventional catheter or the coronary venous leads.
- The tip section of the guidewire is designed for manual shaping. Identify the flexing plane before shaping. Shape in the same plane as that for flexure.
- Guidewire advancement, withdrawal, and torquing should only be carried out when monitored by high quality fluoroscopy
- Do not withdraw or manipulate the coated wire in a metal cannula or sharp-edged object.

ADVERSE EVENTS

Potential adverse events which may result from use of the device include but are not limited to:

- Additional intervention (endovascular or surgical)
- Advancement into the subintimal space
- Allergic reaction
- Antigenic effect
- Arrhythmia
- Cardiac tamponade
- Death
- Dissection
- Embolus
- Emergency Surgery
- Failure to cross occlusion / lesion
- Guidewire entangled/embedded within old stent
- Guidewire lodged in catheter
- Guidewire fracture/separation/peeling
- Hematoma and other access site complications
- Hemorrhage (bleeding)
- Infection
- Irritation to vessel causing vessel spasm
- Myocardial infarction
- No true lumen re-entry
- Perforation

- Pericardial effusion
- Reaction to contrast media
- Thrombus formation
- Vessel dissection or perforation
- Vessel occlusion
- Vessel trauma/injury leading to fistula formation

DIRECTIONS FOR USE

1. Using a sterile, aseptic technique remove the guidewire with the protection spiral from the package and place onto a sterile field.
2. To enhance the gliding performances of the guidewire, flush the protection spiral with saline. A Luer lock connection is attached to the protection spiral.
3. Push gently the proximal end of the guidewire starting from the window of the protection spiral, so that the tip and a portion of the shaft is being exposed out of the Luer lock. Push until the entire proximal end has disappeared into the protection spiral. Pull gently at the distal guidewire shaft, which is then protruding out of the Luer lock until the entire wire is removed out of the protection spiral. Do not grasp the distal tip of the guidewire while removing it from the protection spiral.

The guidewire can be loaded into the interventional catheter system by using one of the two methods described below.

A. Over-the-wire Interventional Catheters

1. Prepare the interventional catheter according to manufacturer's directions.
2. If desired, shape the guidewire tip using standard technique. Do not use a shaping instrument with a sharp edge.
3. Moisten the guidewire with sterile saline to increase the surface lubricity.
4. Place the guidewire into the interventional catheter and advance it until the tip appears just proximal to the distal tip of the interventional catheter.
5. Introduce the interventional catheter/guidewire system through the hemostatic valve on the guiding catheter and advance the system until it is proximal to the tip of the guiding catheter.
6. Secure the interventional catheter within the hemostatic valve being careful not to over-tighten the compression fitting.
7. A torque application device may be applied to the proximal end of the guidewire.
8. Advance the guidewire into the target vessel and across the lesion using fluoroscopy to facilitate proper guidewire placement.
9. Holding the guidewire in position, advance the interventional catheter over the guidewire and into the target lesion.
10. Complete the procedure and remove the guidewire and interventional catheter according to procedural protocol.

B. Rapid Exchange Interventional Catheters

1. Prepare the interventional catheter according to manufacturer's directions.
2. If desired, shape the guidewire tip using standard technique. Do not use a shaping instrument with a sharp edge.
3. Moisten the guidewire with sterile saline to increase surface lubricity.
4. Introduce the guidewire through the hemostatic valve on the guiding catheter using a guidewire introducer.
5. Remove the guidewire introducer by sliding it over the proximal end of the guidewire. Secure the guidewire within the hemostatic valve being careful not to over-tighten the compression fitting.
6. Advance the guidewire into the target vessel and across

the lesion using fluoroscopy to facilitate proper guidewire placement. A torque application device may be applied to the proximal end of the guidewire.

7. Carefully insert the proximal end of the guidewire into rapid exchange interventional catheter and advance it into position according to manufacturer's instructions.
8. Holding the guidewire in position, advance the interventional catheter over the guidewire and into the target lesion.
9. Complete the procedure and remove the guidewire and interventional catheter according to procedural protocol.

The guidewire can be used for positioning compatible over-the-wire (OTW) coronary venous pacing leads in the coronary sinus by using one of the two methods described below.

A. Over-the-wire coronary venous leads

Preload Technique

1. Prepare the lead according to manufacturer's directions.
2. If desired, shape the guidewire tip using standard technique. Do not use a shaping instrument with a sharp edge.
3. Moisten the guidewire with sterile saline to increase surface lubricity.
4. Carefully insert the distal end of the guidewire through the connector pin of the lead.
5. Advance the guidewire until its tip is just proximal to the lead tip.
6. Guide the lead (with guidewire inserted) through the hemostatic valve into the guiding catheter and advance both nearly up to the guiding catheter tip.
7. Tighten the hemostatic valve around the lead to the extent that the lead can just barely be maneuvered.
8. A torque application device may be applied to the proximal end of the guidewire.
9. Under fluoroscopy, advance the guidewire to its target position in the coronary venous system.
10. Maintain guidewire tip position while advancing the lead over the guidewire to the desired position.
11. Carefully remove the guidewire while observing guidewire movement under fluoroscopy.

B. Over-the-wire coronary venous leads

Bare Wire Technique

1. Prepare the lead according to manufacturer's directions.
2. If desired, shape the guidewire tip using standard technique. Do not use a shaping instrument with a sharp edge.
3. Moisten the guidewire with sterile saline to increase surface lubricity.
4. After placing the guide catheter, carefully insert the guidewire through the hemostatic valve that is attached to the guiding catheter using a guidewire introducer.
5. Remove the guidewire introducer by sliding it over the proximal end of the guidewire.
6. Tighten the hemostatic valve around the guidewire to the extent that the wire can just barely be maneuvered.
7. A torque application device may be applied to the proximal end of the guidewire. The torque device must be removed before guiding the lead over the guidewire.
8. Under fluoroscopy, advance the guidewire to its target position in the coronary venous system.
9. Maintain guidewire tip position while advancing the lead over the guidewire to the desired position.
10. Carefully remove the guidewire while observing guidewire movement under fluoroscopy.

CUSTOMER SERVICE INFORMATION

If you have any problems with or questions about Baylis Medical Equipment, contact our technical support personnel.

SYMBOL LEGEND

Explanation of symbols used on the package labels:

	Single use only
	Do not re-sterilize
	Consult Instructions for Use
	STERILE: Sterilized with ethylene oxide gas
	Use by
	Manufacturer
	Date of Manufacture
	Lot Number
	Reference Number
	Keep away from direct sunlight
	Keep dry
	Do not use if package is damaged

LIMITED WARRANTY – DISPOSABLES AND ACCESSORIES

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date. This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

DISCLAIMER AND LIMITATION OF LIABILITY

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY

CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty. Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products - The shelf life of the product

Accessory Products - 90 days from the shipment date

 **Manufactured for:**
BAYLIS MEDICAL COMPANY INC.
5959 Trans-Canada Highway,
Montreal, QC, Canada H4T 1A1
Tel: (514) 488-9801
Worldwide: (800) 850-9801
Web: www.baylismedical.com

OEM:
BRIVANT LTD.
Parkmore West Business Park,
Galway, County Galway
Ireland.
Tel: +353 91 385037
Fax: +353 91 766598

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