



# Initial Experience Using the Radiofrequency Needle Visualization on the Electroanatomical Mapping System for Transseptal Puncture

## INTRODUCTION

- ▶ This series of 42 retrospective consecutive cases evaluates the safety and effectiveness of transseptal puncture (TSP) using a radiofrequency (RF) needle in left-sided ablations with low or no fluoroscopy.

## METHODS

### Visualization setup


- ▶ Pre-procedural contrast enhanced computed tomography or cardiac magnetic resonance were used to create 3D reconstructions of cardiac chambers and vessels.
- ▶ Electroanatomical Mapping (EAM) was performed using CARTO®3 (Biosense Webster) or Rhythmia™ (Boston Scientific) systems and merged with 3D reconstructions using the POLARIS software (Biosense Webster).
- ▶ Transesophageal echocardiography and/or remote magnetic navigation (Stereotaxis Inc) were used in challenging and congenital heart disease (CHD) cases.
- ▶ 3D map of the right atrium (RA) and coronary sinus were acquired using fast anatomical mapping using NaviStar® ThermoCool® Catheter (Bioscience Webster).
- ▶ **NRG™** Transseptal Needle (Baylis Medical\*) was visualized on the EAM map using the **DuoMode™** Extension Cable† (Baylis Medical\*) using the following configuration:
  - A jumper cable (stackable, 2 mm pin) is plugged into ports 1 and 2 on the pin block
  - **DuoMode™** Cable is plugged into the jumper cable in port 1
  - The RF needle was defined as a 2F bipolar catheter, with 2 mm spacing centre-to-centre and 1 mm electrode width/length on the EAM

### Transseptal puncture (TSP)

- ▶ Single or double TSP was performed and 1 or 2 sheaths (**TorFlex™** Transseptal Guiding Sheath, Baylis Medical\* or SL1™ sheath, Abbott) were placed in the left atrium.
  - n=37 cases; first attempt to TSP was made with the **NRG™** Needle
  - n=5 cases; **NRG™** Needle was used after initial attempt with Brockenbrough needle (BN) failed

## RESULTS

- ▶ TSP was achieved 100% successfully with no immediate procedural complications (See Figure 1).



**Figure 1.** The **NRG™** Needle enabled successful TSP with low or no fluoroscopy use in both normal and complex cases with CHD. Number of fluoroleless cases increased with physician experience in **NRG™** Needle visualization on EAM. (Adapted from Guarguagli et al.).

## DISCUSSION AND CONCLUSIONS

- ▶ The prevalence of redo ablations (i.e. fibrotic septa) and adult patients with CHD presents an increasing challenge in TSP.
- ▶ This study demonstrates successful TSP using the **NRG™** Needle in patients with complex and normal interatrial septum anatomies.
- ▶ **NRG™** Needle can be visualized in real time using 3D EAM to reduce/eliminate the need for fluoroscopy.
- ▶ Time saving from more effective TSP using the RF needle offsets the additional time required to map the RA.

\* A wholly-owned subsidiary of Boston Scientific Corporation.

† Consult your mapping system's user manual for connectivity and configuration instructions prior to **DuoMode™** Cable use.

All trademarks are property of their respective owners. Patents Pending and/or issued. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, Contraindications, Warnings, and Instructions For Use can be found in the product labelling supplied with each device or at [www.baylismedical.com](http://www.baylismedical.com).

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.



**Boston  
Scientific**

Advancing science for life™

© 2023 Boston Scientific Corporation  
or its affiliates. All rights reserved.

EP-1562208-AA