
# the FAST TRACK to your therapy delivery in a single solution



Eliminate exchanges



Access with precision




Secure effortless delivery



## **VersaCross**<sup>™</sup> RF Transseptal Solution

Full platform of tools to personalize your solution

### Personalize your solution



#### Brief Summary | **VersaCross**™ RF Wi

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warrings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The Versal roses "RE Wire is indicated for creation of an atrial sental defect in the heart

CONTRAINDICATIONS: The VersaCross MRF Wire is not recommended for use with any conditions that do not require the creation of an atrial septal defect. The Connector Cable is not recommended for use with any other Baylis RF Generator or any other device

WARNINGS: • Laboratory staff and patients can undergo significant x-ray exposure during RF puncture procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. • The VersaCross\*\* RF Wire and Connector Cable provided. Attempts to use it with other connector cables can result in patient complications. • The VersaCross\*\* RF wire and connector Cable provided. Attempts to use it with other connector cables can result in electrocution of the patient and/or operator. • Do not use the VersaCross\*\* RF Wire with electrocautery or electrosurgery generators, connector cables or accessories as attempted use can result in patient and/or operator. • The VersaCross\*\* RF Wire must be used with 0.035\* compatible transseptal sheath and/or dilator devices. Use of incompatible accessory devices may damage the integrity of the VersaCross\*\* RF Wire in sort be used with 0.035\* compatible transseptal puncture use through VersaCross\*\* RF Wire is not intended for use with neonatal patients (i.e. less than one nor fage). Do not attempt to treat neonatal patients with the VersaCross\*\* RF Wire.

PRECAUTIONS: • In order to prevent the risk of ignition, ensure that flammable materials are not present in the room during RF power application. • Careful manipulation of the VersaCross™ RF Wire must be performed to avoid vessel trauma. If resistance is encountered, DO NOT use excessive force to advance or without with VersaCross™ RF Wire or ancillary sheath and/or dilator assembly. Excessive force may lead to bending or kinking of the device limiting advancement and retaction of sheath and/or dilator device. • The Baylis RF Generator is capable of delivering significant electrical pow Patient or operator injury can result from improper handling of the VersaCross™ RF Wire and/or Oil Pelectrode, particularly when operating the device. • During power delivery, the patient should not be allowed to come in contact with ground metal surfaces. • If using electroanatomical mapping guidance, it is recommended to use it along with alternative imaging modality in the event there is loss of visibility of the device.

ADVERSE EVENTS: Adverse events that may occur while creating an atrial septal defect include: • Tamponade • Sepsis/Infection • Thromboembolic episodes • Vessel perforation • Atrial Fibrillation • Myocardial Infarction • Vessel spasm • Sustained arrhythmias • Atrial Flutter • Hemorrhage • Vascu thromboss • Perforation of the myocardium • Hematoma • Additional Surgical Procedure • Wire entrapment/ entanglement • Foreign body/ wire fracture
EP-1504711-AA

#### Brief Summary | **VersaCross**™ Transseptal Dilator

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The VersaCross® Transsental Dilator is indicated for introducing various cardiovascular catheters to the heart including the left side of the heart through the interatrial sentum.

CONTRAINDICATIONS: There are no known contraindications for this device

WARNINGS: • Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure - The Versacross\*\* Steerable Sheath kit is intended for single patient use only. Do not attempt to sterifize and reuse the Versacross\*\* Steerable Sheath kit. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Sallure to follow this instruction may recurring the patient and in patient to another Sallure to follow this instruction may recurring the patient and in patient to another sallure to follow this instruction may recurring the patient to another sallure to follow this instruction may recurring the patient of the patient the patient may recurrent the patient of the patient may recurrent the patient may recurrent

PRECAUTIONS: • Careful manipulation must be performed to avoid cardiac damage or tamponade. Sheath, dilator and guidewire advancement should be done under fluoroscopic guidance if resistance is encountered. DO NOT use excessive force to advance or withdraw the device

ADVERSE EVENTS: Adverse events that may occur while using the Versacross. Sheath include: • Infection • Air embolus • Local nerve damage • Vasovagal reaction • Dissection • Vessel spasm • AV fistula formation • Atrial septal defect • Pseudoaneurysm • Aortic puncture • Arrhythmias • Perforati and/or tamponade • Hemotronage • Catherorrape • Coronary artery spasm and/or damage • Wasovagal reaction • Pacemaker/defibrillator lead displacement • Pulmonary edema • Coronary artery spasm and/or damage • Vessel tamma • Pericardial/pleural effusion

#### Brief Summary | **VersaCross™** Transseptal Sheat

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warning Precautions, Adverse Events, and Operator's Instructions.

IDICATIONS FOR USE: The VersaCross \*\* Transseptal Sheath kit is used for the percutaneous introduction of various types of cardiovascular catheters and guidewires to all heart chambers, including the left atrium via transseptal erforation / puncture.

#### NTRAINDICATIONS: There are no known contraindications for this device.

ARANINGS: • Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well a receased risk for somatic and generic effects. Therefore, adequate measures must be taken to minimize this exposure. The use of enchardingaphly is recommended. • The Versacross\*\* Transsepatl Sheath kit is intended for single attent use only. Do not attempt to sterilize and reuse the Versacross\*\* Transsepatal Sheath kit. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another. • Do not attempt direct encurations of the sheath without the dilator as this may cause vessel injury. • Careful manipulation must be performed to avoid cardiac damage or tamponade. Sheath advancement should be done under fluoroscopic undance. Education graphic guidatence is also recommended.

AUTIONS: • Careful manipulation must be performed to avoid cardiac damage or tamponade. Sheath, dilator and guidewire advancement should be done under fluoroscopic guidance. If resistance is encountered, DO NOT use sive force to advance or withdraw the device. • The VersaCross™ Transseptal Sheath and Dilator are compatible with transseptal devices juidewires. 035' or smaller. • The VersaCross™ Transseptal Sheath and Dilator are compatible with transseptal needles such as the "NRG™ Transseptal Needle".

VERSE EVENTS: Adverse events that may occur while using the VersaCross<sup>™</sup> Transseptal Sheath kit include: • Infection • Air embolus • Local nerve damage • Hemorrhage • Embolic events • Vessel spasm • AV fistula formation • also pela defect • Pseudoaneurysm • Perforation and/or tamponade • Arrhythmilas • Pericardial/pleural effusion • Hematoma • Vessel trauma • Valve damage • Catheter entrapment 6 of the Company of the Comp

Boston Scientific is a Global Company. Please note that model numbers, indications, contraindications, warnings and specifications may differ depending on geographic region. Not all information displayed in this brochure may be licensed in accordance with Canadian law. Please contact your Boston Scientific representative for local labeling, product specifications and licensed model numbers. CAUTION: The law restricts this device to sale by or on the order of a physician. Rx only. Indications, Contraindications, Warnings, and Instructions For Use can be found in the product labelling supplied with each device or at www.bylismedical.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries.

# Scientific Advancing science for life

Baylis Medical Company Inc. 5959 Trans-Canada Highway Montreal. OC Canada H4T 1A1

www.baylismedical.com


nfo@baylismedical.com

eral Inquiries ) 488-9801

© 2023 Boston Scientific Corporation or its affiliates. All rights reserved. EP-1545803-AA







# TRANSSEPTAL REIMAGINED

Fast track access-to-delivery with a single solution

# the only EXCHANGELESS\* solution for access-to-delivery of left heart therapy devices


# VERSATILITY

### **ELIMINATE EXCHANGES**

Left heart access using a SINGLE SOLUTION from start to finish









Rewire and drop down to optimize transseptal location WITHOUT THE HASSLE of exchanging a needle

ersaCross™ RF Wire can be used, without exchanges, as a guidewire, as a ransseptal puncture device or as an exchange rail for delivering therapy sheath

### ACCESS WITH PRECISION

PRECISION RF PUNCTURE TECHNOLOGY to optimize transseptal location for any anatomy





## TRUform<sup>™</sup> SHAPEABLE TECHNOLOGY

with true-to-form curve retention





Snap-fit to maintain handsfree sheath-dilator

### SECURE EFFORTLESS DELIVERY

Instantly gain and maintain access WITHOUT EXCHANGES



Sturdy 0.035" rail to deliver therapy devices with confidence


Flexible tip cushions against left atrial wall and reduces risk of perforation

Know where you are at all times with

**OMNIviz**<sup>™</sup> Technology



Visualize your entire solution on fluoroscopy



Reliably locate your devices on ultrasound to reduce reliance on fluoroscopy



Visibly confirm position of RF tip within dilator



Track and mark RF tip position on your mapping system

VersaCross™ RF Transseptal Solution