## Instructions for Use DuoMode™ Cable

English.....	1
Français.....	2
Italiano.....	3
Deutsch.....	4
Dutch.....	6
Espanol.....	7
Portugues.....	8
Čeština.....	10
Dansk.....	11
Suomi.....	12
Norsk.....	13
Svenska.....	14
Slovenčina.....	15
Română.....	16
Hrvatski.....	18
Magyar.....	19
Türkçe.....	20
Ελληνικά.....	21
Русский.....	22

Baylis Medical Company Inc.  
5959 Trans-Canada Highway  
Montreal, Quebec, Canada, H4T 1A1  
Tel: (514) 488-9801/ (800) 850-9801 Fax: (514) 488-7209  
www.baylismedical.com

EU Authorized Representative:  
Quality First International OU  
Laki 30, 12915 Tallinn  
Estonia  
Telephone: +372 610 41 96  
Email: [enquiries@qualityfirstint.ee](mailto:enquiries@qualityfirstint.ee)

© Copyright Baylis Medical Company Inc., 2023

Baylis Medical and the Baylis Medical logo are trademarks of Baylis Medical Technologies Inc.

**English**  
Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

**Caution:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.  
**Notice (EU):** Any serious incident as defined in the EU Medical Device Regulation (2017/745) that has occurred in relation to the DuoMode Cable should be reported to Baylis Medical Company Inc. and the competent authority of the EU Member State in which the physician and/or patient is established.

**I. DEVICE DESCRIPTION**  
The reusable DuoMode Cable is an extension cable that is used with Baylis Medical approved radiofrequency puncture devices, the Baylis Medical Company Radiofrequency Puncture Generator (RFP-100A) and diagnostic equipment. Detailed information concerning the Generator is contained in a separate manual that accompanies the Generator (RFP-100A Generator Instructions for Use). In addition, detailed information concerning the RF puncture devices is contained in separate manuals that accompany these devices.

The DuoMode Cable has a four-pin connector on a cable that mates with the RFP Generator and a puncture device connector port that accepts other Baylis RFP connector cables that facilitates connection to the puncture device. The DuoMode connector cable also has a Diagnostic Equipment Connector that consists of a protected 2mm pin.

**II. INDICATIONS FOR USE**  
The DuoMode Cable is intended to serve as an extension cable that is used with the Baylis Medical radiofrequency puncture devices, the Baylis Medical Company Radiofrequency Puncture Generator and diagnostic equipment.

**III. CONTRAINDICATIONS**  
The DuoMode Cable is not recommended for use with any other RF generator.

**IV. WARNINGS**

- The DuoMode Cable is a reusable device. Failure to properly clean the device can cause patient injury and/or the communication of infectious disease(s) from one patient to another.
- The DuoMode Cable must only be used with Baylis RF Puncture Generators and RF puncture devices. Attempts to use it with other RF Generators and devices can result in electrocution of the patient and/or operator.
- Laboratory staff and patients can undergo significant x-ray exposure during radiofrequency puncture procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure.

**V. PRECAUTIONS**

- Do not attempt to use the DuoMode Cable or ancillary equipment before thoroughly reading the accompanying Instructions for Use.
- Puncture procedures should be performed only by physicians thoroughly trained in the techniques of RF powered puncture and use of diagnostic medical equipment in a fully equipped catheterization laboratory.
- Visually inspect the DuoMode Cable to ensure there is no degradation, such as cracking or damage to the insulating material. Do not use the cable if there is any damage.
- The DuoMode Cable is intended for use with RF puncture devices only.
- Never disconnect the DuoMode Cable from the RF Puncture Generator while the Generator is delivering RF power.
- Never disconnect the DuoMode Cable from the Generator by pulling on the cable. Failure to disconnect the cable properly may result in damage to the cable.
- Do not twist the DuoMode Cable while inserting or removing it from the Isolated Patient Connector on the Generator. Twisting the cable may result in damage to the pin connectors.
- Do not bend the cable. Excessive bending or kinking of the cable may damage the integrity of the cable and may cause patient injury. Care must be taken when handling the cable.

- Take precautions to limit the effects that the electromagnetic interference (EMI) produced by the Generator may have on the performance of other equipment. Check the compatibility and safety of combinations of other physiological monitoring and electrical apparatus to be used on the patient in addition to the Generator.
- Adequate filtering must be used to allow continuous monitoring of the surface electrocardiogram (ECG) during RF power applications.
- During power delivery, the patient should not be allowed to come in contact with grounded metal surfaces.
- In order to prevent the risk of ignition make sure that flammable material is not present in the room during RF power application.

Baylis Medical Company relies on the physician to determine, assess and communicate to each individual patient all foreseeable risks of the Baylis Medical Radiofrequency Puncture System.

**VI. ADVERSE EVENTS**  
The RFP-100A Generator Instructions for Use should be consulted for any adverse events.

**VII. PRODUCT SPECIFICATIONS**


Model Number	RFX-BAY-DUO-100A
Generator Connector	4-pin (Plug)
Generator Connector Cable Colour	Black
Diagnostic Equipment Connector	Touchproof 2mm Connector (DIN 42802)
Diagnostic Equipment Cable Colour	Red
Puncture Device Connector Port	4-pin (receptacle)

**VIII. INSPECTION PRIOR TO USE**  
Perform the following checks before the patient is presented for the procedure. These tests will allow you to verify that the equipment you will use is in proper working order. Do these tests in a sterile environment. Do not use defective equipment.

KEY ITEMS	QUESTION?	WARNINGS AND EXPLANATIONS
Visual Check	Have you done a visual check on the entire system?	Ensure the connectors, the enclosure and the cable have no visible damage, such as discoloration, cracks, label fading, cable splice, or kinks. Do not use damaged equipment.

**IX. EQUIPMENT REQUIRED**  
Puncture procedures should be performed in a specialized clinical setting which may be equipped with a fluoroscopy unit, radiographic table, physiologic recorder, emergency equipment and instrumentation for gaining vascular access.

**X. DIRECTIONS FOR USE**  
Once the RF puncture device is properly positioned at the puncture site, and the Generator is properly set up (following the instructions in the Generator Instructions for Use), the DuoMode Cable can be used to connect the catheter or wire to the Generator.

**Connecting the DuoMode cable**

- Connect the Generator connector end (black cable) of the DuoMode Cable to the isolated patient connector port on the RF Puncture Generator as per the Generator Instructions for Use. The DuoMode Cable Generator Connector uses a circular connector, keyed for proper alignment. Gently line up the key features of the connector and the connector pins with the socket and push in until the connector fits firmly into the socket. Any attempt to connect the cable otherwise will damage the pins on the connector.
- Connect the Diagnostic Equipment Connector (red cable) to the input of the diagnostic equipment.
- Connect the Generator connector end of the RFP Connector Cable (used to connect the puncture device) to the Puncture Device Connector Port on the DuoMode cable. Follow the RFP Connector Cable Instructions for Use.

**Note:** The generator connector of the RFP Connector cable can plug directly into the isolated patient connector on the Generator, if the DuoMode Cable is not used.

- Do not use excessive force in connecting any of the cables. Use of excessive force may result in damage to the connector pins.
- Connect the device connector end of the RFP Connector Cable to the RF Puncture Device, according to the RFP Connector Cable Instructions for Use.

- Switching the DuoMode Cable**
- To connect the puncture device to the Diagnostic Equipment, press the rocker switch down towards the Mapping symbol.
  - To connect the puncture device to the RF Puncture Generator, press the rocker switch down towards the Generator symbol. Do not activate RF delivery on the generator while the switch is set to the Mapping setting. Do not change the rocker switch position while RF is being delivered.

**Disconnecting the DuoMode Cable**

- To disconnect the puncture device from the RFP Connector Cable, follow the RFP Connector Cable Instructions for Use.
- To disconnect the RFP Connector Cable from the DuoMode cable, follow the RFP Connector Cable Instructions for Use.
- To disconnect the Generator Connector from the generator, while firmly holding the Generator Connector housing, gently pull the Generator Connector straight out of the isolated patient connector socket on the Generator.
- To disconnect the Diagnostic Equipment, grasp the Diagnostic Equipment Connector firmly and gently pull it straight out of the socket.

**XI. CLEANING AND STERILIZATION INSTRUCTIONS**  
The DuoMode cable is a non-body contact device and therefore is NOT sterile and CANNOT be sterilized. If cleaning is necessary, the surface of the DuoMode cable can be cleaned with a damp low lint cloth with non-abrasive detergent dissolved in water. Dry the surface after wiping down. Do not spray or pour liquids directly on the DuoMode Cable.

**XII. DISPOSAL OF WASTE**  
The DuoMode Cable is reusable. However, when the DuoMode Cable is determined to be unusable, dispose of the DuoMode Cable in compliance with standard hospital procedures.

**XIII. CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION**  
If you have any problems with or questions about Baylis Medical Equipment contact our technical support personnel.

- NOTES:**
- In order to return products you must have a return authorization number before shipping the products back to Baylis Medical Company.
  - Baylis Medical will not accept any piece of used equipment without a sterilization certificate. Ensure that any product being returned to Baylis Medical has been cleaned, decontaminated and sterilized as per user instructions before returning it for warranted service.

**XIV. TROUBLESHOOTING**  
The following table is provided to assist the user in diagnosing potential problems.

PROBLEM	COMMENTS	TROUBLESHOOTING
<b>Generator Alert Messages</b>	In order to successfully puncture tissue using radiofrequency energy, the entire system must be connected and all devices must be in good working order.	Ensure that all connections are made: - puncture device to connector cable - connector cable to DuoMode cable - DuoMode cable to generator - generator to power outlet - generator to grounding pad Visually inspect the catheter/wire or cable for damage. Immediately discard any damaged equipment. If the problem persists discontinue use. For error/alert messages encountered while attempting puncture, refer to the operator's manual that accompanies the Generator. If errors persist, attach a new connector cable. If this solves the problem, discard the damaged connector cable.
<b>DuoMode Cable does not fit into the Isolated Patient Connector on the front panel of the generator</b>	The connectors are designed to connect in a specific way for safety reasons. If the connector "keys" are out of line, the connectors won't fit together	Check that the connector keys are lined up in the proper orientation. Ensure that the connectors are clean and unobstructed.

#### XV. LABELING AND SYMBOLS

	Manufacturer		Caution
	Consult Instructions for Use		Lot Number
<b>REF</b>	Model number		Keep Away From Sunlight
	Date of Manufacture		Do Not Use if Packaging is Damaged
	EU Authorized Representative		Medical Device
	To Generator		To Diagnostic Equipment
<b>Rx ONLY</b>	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.		Cable
<b>USE WITH: BMC RFP-100A GENERATORS ONLY</b>		Use with: BMC RFP-100A generators only	
<b>DO NOT STEAM STERILIZE</b>		Do Not Steam Sterilize	
	Only for EU member states: Use of this symbol indicates that the product must be disposed of in a way that complies with local and national regulations. For questions regarding recycling of this device please contact your distributor		

#### XVI. LIMITED WARRANTY – Disposable and Accessories

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date.

This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

#### DISCLAIMER AND LIMITATION OF LIABILITY

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HERUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products	The shelf life of the product
Accessory Products	90 days from the shipment date

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