Instructions for Use

VersaCross Connect™ Transseptal Dilator

English.....	1
Français	3
Deutsch	5
Nederlands	7
Italiano	9
Español.....	11
Português	13
Čeština.....	15
Dansk	17
Suomi.....	19
Norsk	21
Svenska	23

Baylis Medical Company Inc.
 5959 Trans-Canada Highway
 Montreal, Quebec, Canada, H4T 1A1
 Tel: (514) 488-9801/ (800) 850-9801
 Fax: (514) 488-7209
 www.baylismedical.com

©2023 Baylis Medical Company Inc. or its affiliates.

English

All trademarks are property of their respective owners.

Carefully read all instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications. Baylis Medical Company relies on the physician to determine, assess, and communicate to each individual patient all foreseeable risks of the procedure. CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

DEVICE DESCRIPTION

The VersaCross Connect Transseptal Dilator is designed for safe and easy catheterization and angiography of specific heart chambers and locations. The dilator provides torque control and is flexible. The dilator features a tapered tip and a shaft that can be reshaped manually. The echogenic shaft and tip and radiopaque tip maximize visualization of the dilator during manipulation.

The VersaCross Connect Transseptal Dilator is for use with a 13F ID FARADRIVE™ Steerable Sheath which is 74cm in length, specifically, models: M004PF21M402 (USA), M004PFCE21M402 (Canada).

Carefully read the applicable FARADRIVE Steerable Sheath instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

INDICATIONS FOR USE

The VersaCross Connect Transseptal Dilator is used for the percutaneous introduction of various types of cardiovascular catheters and guidewires to all heart chambers, including the left atrium via transseptal perforation / puncture.

United States: The VersaCross Connect Transseptal Dilator is indicated for use in procedures where access to the left atrium via the transseptal technique is desired.

WARNINGS

- Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. The use of echocardiography is recommended.

- The VersaCross Connect Transseptal Dilator is intended for single patient use only. Failure to do so may result in patient complications.
- Do not attempt to sterilize and reuse the VersaCross Connect Transseptal Dilator. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another.
- The VersaCross Connect Transseptal Dilator is supplied STERILE using an ethylene oxide process. Do not use if the package is damaged.
- Care should be taken to ensure that all air is removed from the dilator before infusing through the proximal hub.
- Care should be taken when inserting or removing the dilator from access sheaths or introducer sheaths.
- Manual shaping of the distal curve shall be done with smooth motions along the curve. Do not use excessive force and/or pressure when reshaping.
- Care should be taken when inserting or removing compatible guidewires from the dilator lumen.
- Do not attempt direct percutaneous insertion of the dilator without a guidewire as this may cause vessel injury.
- Careful manipulation must be performed to avoid cardiac damage or tamponade. Dilator advancement should be performed under imaging guidance. Fluoroscopic or echocardiographic imaging is recommended. If resistance is encountered, DO NOT use excessive force to advance or withdraw the device.

PRECAUTIONS

- Do not attempt to use the VersaCross Connect Transseptal Dilator before thoroughly reading the accompanying Instructions for Use.
- The sterile barrier system and dilator should be visually inspected prior to use. Do not use if the sterile barrier integrity or device have been compromised or damaged.
- Only physicians thoroughly trained in the techniques of the approach to be used should perform interventional procedures.
- The VersaCross Connect Transseptal Dilator is compatible with introducer sheaths 13F or larger.
- The VersaCross Connect Transseptal Dilator is for use with the 13F ID FARADRIVE Steerable Sheath, which is 74cm in length.
- The VersaCross Connect Transseptal Dilator is compatible with 0.035" transseptal devices and guidewires or smaller.
- The VersaCross Connect Transseptal Dilator is NOT compatible with transseptal needles such as the "NRG™ Transseptal Needle".

CONTRAINDICATIONS

There are no known contraindications for this device.

SPECIAL STORAGE AND/OR HANDLING INSTRUCTIONS

Keep away from sunlight.

ADVERSE EVENTS

Adverse events that may occur while using the VersaCross Connect Transseptal Dilator include:

Infection	Air embolus
Local nerve damage	Vessel trauma
Vessel spasm	Pseudoaneurysm
AV fistula formation	Atrial septal defect
Arrhythmias	Perforation and/or tamponade
Hematoma	Hemorrhage
Catheter entrapment	Embolic events
Valve damage	Pericardial/pleural effusion

PREPARATION FOR USE

Prior to use of the VersaCross Connect Transseptal Dilator, the device should be carefully examined for damage or defects, as should all equipment used in the procedure. Do not use defective equipment. Do not reuse the device.

SUGGESTED DIRECTIONS FOR USE

- Carefully read all instructions prior to use. Failure to do so may result in complications.
- Carefully read the applicable FARADRIVE Steerable Sheath instructions prior to use. Failure to do so may result in complications.
- Carefully read the applicable compatible guidewire instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications.
- Use of VersaCross Connect is expected to reduce the number of exchanges in the procedure resulting in a more efficient transseptal puncture. This should be taken into account when estimating the timing for heparin administration to ensure appropriate ACT levels after transseptal puncture.
- Dilator may be used with a compatible sheath or as a standalone device to facilitate access to the left atrium via transseptal puncture.
- Thoroughly flush the dilator with heparinized saline solution prior to use.
- Gain access to the right femoral vein using standard methods.
- Care should be taken when inserting or removing the dilator from the venous cutaneous puncture site; note that a compatible access introducer sheath may be used if desired. Refer to the compatible introducer sheath's Instructions for Use for details and directions.
- Introduce a compatible guidewire through the vasculature access point and advance to required depth.

- Dilator can be inserted fully into the compatible access sheath (if sheath is used) and a manual curve may be added to the dilator or dilator and sheath assembly prior to insertion into the body.
- Thread the dilator and compatible access sheath (if sheath is used) over the guidewire, allowing the device to twist freely while under imaging guidance such as fluoroscopy or echocardiography up to the SVC. If resistance is encountered, DO NOT use excessive force to advance or withdraw the dilator over the guidewire. Determine the cause of resistance before proceeding.
- Use standard technique to position the dilator/guidewire assembly or the sheath/dilator/guidewire assembly (if sheath is used) into the desired heart chamber.
- If transseptal puncture is required, refer to the Instructions for Use of the transseptal puncture device.
- Ensure the dilator is clear of air. To aspirate blood, use the dilator hub.
- Monitor the location of the radiopaque tip frequently under imaging guidance, such as fluoroscopy or echocardiography.
- Deliver a continuous heparinized solution infusion or aspirate periodically. This may help reduce the risk of thromboembolic complications due to thrombus formation, as there may be a possibility of thrombus development at the distal dilator tip or inside the dilator lumen. Also aspirate when removing the transseptal device or dilator.
- After removal of the dilator, use standard technique to achieve hemostasis.

CLEANING AND STERILIZATION INSTRUCTIONS

Do not clean or re-sterilize the VersaCross Connect Transseptal Dilator. The VersaCross Connect Transseptal Dilator is intended for single use only.

DISPOSAL OF WASTE

Treat the used device(s) as biohazardous waste and dispose of in compliance with standard hospital procedures.

CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION


If you have any problems with or questions about Baylis Medical Equipment, contact our technical support personnel.

Baylis Medical Company Inc.
 5959 Trans-Canada Highway
 Montreal, Quebec, Canada, H4T 1A1
 Phone: (514) 488-9801 or (800) 850-9801
 Fax: (514) 488-7209
 www.baylismedical.com

NOTES:

1. In order to return products, you must have a return authorization number before shipping the products back to Baylis Medical Company. Product Return Instructions will be provided to you at this time.
2. Ensure that any product being returned to Baylis Medical has been cleaned, decontaminated and/or sterilized as indicated in the Product Return Instruction before returning it for warrantied service. Baylis Medical will not accept any piece of used equipment that has not been properly cleaned or decontaminated as per the Product Return Instructions.

LABELING AND SYMBOLS

	Manufacturer	Rx ONLY	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
	Sterile using ethylene oxide		Single Use – Do not reuse
	Use By		Batch Code
	Caution		Do not use if package is damaged and consult instructions for use.
	Consult Instructions for Use		Keep Away from Sunlight
	Model number		Do Not Re-sterilize
	Non-Pyrogenic		Maximum guidewire outside diameter that can be used with this device

LIMITED WARRANTIES

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to

BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date. This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

DISCLAIMER AND EXCLUSION OF OTHER WARRANTIES

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

LIMITATION OF LIABILITY FOR DAMAGES

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products	The shelf life of the product
Accessory Products	90 days from the shipment date

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