Instructions for Use

ExpanSure™ Large Access Transseptal Dilator

[English	1]
[Français	3]
[Deutsch	4]
[Nederlands	6]
[Italiano	8]
[Español	9]
[Português	11]
[Čeština	13]
[Dansk	14]
[Suomeksi	16]
[Norsk	17]
[Svenska	19]

Baylis Medical Company Inc.
5959 Trans-Canada Highway
Montreal, Quebec, Canada, H4T 1A1
Tel: (514) 488-9801/ (800) 850-9801
Fax: (514) 488-7209
www.baylismedical.com

EU Authorized Representative:
Quality First International OÜ
Laki 30, 12915 Tallinn
Estonia
Telephone: +372 610 41 96
Email: enquiries@qualityfirstint.ee

© Copyright Baylis Medical Company Inc., 2019 - 2021

The Baylis Medical and the Baylis Medical logo are trademarks of Baylis Medical Technologies Inc.

English

Carefully read all instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications. Baylis Medical Company relies on the physician to determine, assess and communicate to each individual patient all foreseeable risks of the procedure.

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

DEVICE DESCRIPTION

The ExpanSure™ Large Access Transseptal Dilator consists of two components: a dilator and a J-tipped Mechanical guidewire.

The ExpanSure Large Access Transseptal Dilator is designed for safe and easy catheterization and angiography of specific heart chambers and locations. The dilator provides superior torque control and is flexible. The radiopaque and echogenic tip maximizes visualization of the dilator during manipulation.

The J-tipped Mechanical Guidewire, hereafter referred to as the "guidewire", is comprised of a stainless-steel core with a flexible, spiral shaped PTFE coated steel coil along the full length of this device. The dilator shaft and guidewire are in their entirety coated with hydrophobic lubricious coatings for smoother device manipulation. No pre-conditioning is required for these coatings.

INDICATIONS FOR USE

The ExpanSure Large Access Transseptal Dilator is used for the percutaneous introduction of various types of cardiovascular catheters and guidewires to all heart chambers, including the left atrium via transseptal perforation / puncture.

United States: The ExpanSure Large Access Transseptal Dilator is indicated for use in procedures where access to the left atrium via the transseptal technique is desired.

WARNINGS

- Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. The use of echocardiography is recommended.
- The ExpanSure Large Access Transseptal Dilator's shaft is coated with a lubricious coating. The following warnings must be considered:
 - Use with sheaths smaller than the size(s) listed in the section below may result in a tight fit that affects device performance, including coating integrity.
 - Excessive wiping and/or wiping with a dry gauze may damage the coating.
 - Excessive manual bending and/or shaping of the shaft may affect the coating integrity.
- The Mechanical Guidewire is coated with a lubricious coating. The following warnings must be considered:
 - Use with incompatible introducers or dilators may affect device performance and integrity, including coating integrity
 - Excessive manual bending and/or shaping of the device may affect the coating integrity
- The ExpanSure Large Access Transseptal Dilator and accompanying guidewire is intended for single patient use only. Do not attempt to sterilize and reuse the ExpanSure Large Access Transseptal Dilator or the accompanying guidewire. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Failure to do so may result in patient complications.
- The ExpanSure Large Access Transseptal Dilator and accompanying guidewire are supplied STERILE using an ethylene oxide process. Do not use if the package is damaged.
- Care should be taken to ensure that all air is removed from the dilator before infusing through the proximal hub.
- DO NOT attempt to insert or retract the guidewire through a metal cannula or a percutaneous needle, which may damage the guidewire and may cause patient injury
- Careful manipulation must be performed to avoid cardiac damage or tamponade. Dilator and guidewire advancement should be done under imaging guidance such as fluoroscopy or echocardiography.

PRECAUTIONS

- Careful manipulation must be performed to avoid cardiac damage, or tamponade. Dilator and guidewire advancement should be done under imaging guidance such as fluoroscopy or echocardiography. If resistance is encountered, DO NOT use excessive force to advance or withdraw the device.
- The sterile packaging, dilator, and guidewire should be visually inspected prior to use. Do not use the device if it has been compromised or damaged.
- Do not attempt to use the ExpanSure Large Access Transseptal Dilator or the accompanying guidewire before thoroughly reading the accompanying Instructions for Use.
- Only physicians thoroughly trained in the techniques of the approach to be used should perform interventional procedures.
- Note product "Use By" date.

- The “ExpanSure Large Access Transseptal Dilator” is compatible with the “NRG™ Transseptal Needle”.
- The “ExpanSure Large Access Transseptal Dilator” is compatible with 12.5F introducer sheaths.
- The “ExpanSure Large Access Transseptal Dilator” is compatible with .032” and .035” transseptal devices.
- The “ExpanSure Large Access Transseptal Dilator” is compatible with .032” and .035” guidewires

CONTRAINDICATIONS

- There are no known contraindications for this device.

SPECIAL STORAGE AND/OR HANDLING INSTRUCTIONS

- Keep away from sunlight.

ADVERSE EVENTS

Adverse events that may occur while using the ExpanSure Large Access Dilator include:

- | | |
|----------------------|------------------------------|
| Infection | Air embolus |
| Local nerve damage | Vessel trauma |
| Vessel spasm | Pseudoaneurysm |
| AV fistula formation | Atrial septal defect |
| Arrhythmias | Perforation and/or tamponade |
| Hematoma | Hemorrhage |
| Catheter entrapment | Embolic events |
| Valve damage | Pericardial/pleural effusion |

PREPARATION FOR USE

Prior to use of the ExpanSure Large Access Transseptal Dilator, the dilator and guidewire should be carefully examined for damage or defects, as should all equipment used in the procedure. Do not use defective equipment. Do not reuse the device.

Intracardiac puncture procedures should be performed in a sterile environment in a specialized clinical setting equipped with a fluoroscopy unit, radiographic table, echocardiographic imaging, physiologic recorder, emergency equipment and instrumentation for gaining vascular access.

SUGGESTED DIRECTIONS FOR USE

- Carefully read all instructions prior to use. Failure to do so may result in complications.
- Thoroughly flush the guidewire and dilator with heparinized saline solution prior to use.
- Perform a standard vein puncture of the right femoral vein using an access needle (not supplied).
- Introduce the guidewire through the vasculature access point and advance to the required depth. If resistance is encountered, DO NOT use excessive force to advance or withdraw the guidewire. Determine the cause of resistance before proceeding.
- Enlarge the cutaneous puncture site as necessary.
- Note that a compatible access introducer sheath may be used at the venous cutaneous puncture site if desired. Refer to the compatible introducer sheath’s Instructions for Use for details and directions.
- Thread the dilator over the guidewire using a slight twisting motion into the superior vena cava (SVC) under imaging guidance such as fluoroscopy or echocardiography. If resistance is encountered, DO NOT use excessive force to advance or withdraw the dilator/sheath assembly over the guidewire. Determine the cause of resistance before proceeding.
- Use standard technique to position the dilator into the desired heart chamber.
- If transseptal puncture is required, refer to the Instructions for Use of the transseptal puncture device.
- Ensure the dilator is clear of air. To aspirate blood, use the dilator hub.
- Monitor the location of the radiopaque tip frequently under imaging guidance such as fluoroscopy or echocardiography.
- Deliver a continuous heparinized solution infusion or aspirate periodically. This may help reduce the risk of thromboembolic complications due to thrombus formation, as there may be a possibility of thrombus development at the distal dilator tip or inside the dilator lumen. Also aspirate when removing the transseptal perforating device or dilator.
- After removal of the dilator, use standard technique to achieve hemostasis.

CLEANING AND STERILIZATION INSTRUCTIONS

Do not clean or re-sterilize the ExpanSure Large Access Transseptal Dilator or the accompanying guidewire. The ExpanSure Large Access Transseptal Dilator and accompanying guidewire are intended for single use only.

DISPOSAL OF WASTE

Treat the used device(s) as biohazardous waste and dispose of in compliance with standard hospital procedures.

CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION

If you have any problems with or questions about Baylis Medical Equipment contact our technical support personnel.

Baylis Medical Company Inc.
 5959 Trans-Canada Highway
 Montreal, Quebec, Canada, H4T 1A1
 Phone: (514) 488-9801 or (800) 850-9801
 Fax: (514) 488-7209
 www.baylismedical.com

NOTES:

1. In order to return products, you must have a return authorization number before shipping the products back to Baylis Medical Company. Product Return Instructions will be provided to you at this time.
2. Ensure that any product being returned to Baylis Medical has been cleaned, decontaminated and/or sterilized as indicated in the Product Return Instruction before returning it for warranted service. Baylis Medical will not accept any piece of used equipment that has not been properly cleaned or decontaminated as per the Product Return Instructions.

LABELING AND SYMBOLS

	Manufacturer	Rx ONLY	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
	Sterile using ethylene oxide		Do not reuse
	Use-by date		Batch Code
	Non-Pyrogenic		Do Not Re-Sterilize
	Caution		Do Not Use if Packaging is Damaged
	Consult Instructions for Use		Keep Away from Sunlight
	Model number		EU Authorized Rep

LIMITED WARRANTIES

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date. This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC’s Limited Warranty shall not apply to BMC products which have been reesterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC’s instructions.

DISCLAIMER AND LIMITATION OF LIABILITY

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER’S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER’S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR

DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products	The shelf life of the product
Accessory Products	90 days from the shipment date

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