Instructions for Use**ProTrack™ Pigtail Wire**

English.....	1
Français.....	3
Deutsch.....	5
Nederlands.....	8
Italiano.....	10
Español.....	13
Português.....	15
Čeština.....	18
Dansk.....	20
Suomi.....	22
Norsk.....	24
Svenska.....	26

Baylis Medical Company Inc.
5959 Trans-Canada Highway
Montreal, Quebec, Canada, H4T 1A1
Tel: (514) 488-9801/ (800) 850-9801
Fax: (514) 488-7209
www.baylismedical.com

EU Authorized Representative:
Quality First International OÜ
Laki 30, 12915 Tallinn
Estonia
Telephone: +372 610 41 96
Email: enquiries@qualityfirstint.ee

EU Importer:
MedEnvoy
Prinses Margrietplantsoen 33 – Suite 123
2595 AM The Hague
The Netherlands

© Copyright Baylis Medical Company Inc., 2007 - 2022

The Baylis Medical logo, Baylis Medical and ProTrack are trademarks and/or registered trademarks of Baylis Medical Company Inc. in the United States of America and/or other countries.

English

Carefully read all instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

Baylis Medical Company relies on the physician to determine, assess and communicate to each individual patient all foreseeable risks of the procedure.

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

I. DEVICE DESCRIPTION

Baylis Medical *ProTrack™* Pigtail Wires, from hereafter referred to as "pigtail wires", are comprised of a stainless steel core wire with a flexible, spiral shaped stainless steel coil on the distal end. The Pigtail wire is offered in two lengths including, 175 and 230 cm. The product is supplied sterile for single use only.

NOTE: Refer to product labeling for specific model information regarding pigtail wire diameters, lengths and tip shapes.

II. INDICATIONS FOR USE

The *ProTrack™* Pigtail Wires are intended for use in percutaneous transeptal procedures to introduce and position catheters and other interventional devices within the left heart. The device is not intended for use in the coronary arteries.

III. WARNINGS

- The pigtail wire is intended for single use only. DO NOT resterilize and/or reuse.
- The pigtail wire should only be used by a physician trained in percutaneous, intravascular techniques and procedures. It is important to follow the instructions for use prior to using this product.
- The guidewire should only be handled by users who are familiar with aseptic handling techniques
- The pigtail wire is supplied sterile and non-pyrogenic in an unopened package and should be used under sterile conditions. DO NOT use if the package has been damaged or opened.
- DO NOT use the pigtail wire after the expiration date indicated on the label. Discard pigtail wires that exceed the expiration date.
- Confirm ancillary devices are compatible with the pigtail wire outer diameter before use.
- Always exercise care while advancing or withdrawing the pigtail wire to prevent wire breakage, kinking, coil separation or damage. Wire damage may cause complications which may require additional interventions.
- DO NOT push, auger, withdraw or torque a pigtail wire against resistance until the cause of the resistance has been determined. Applying excessive force against unexpected resistance may cause damage to the pigtail wire, interventional device and/or vessel/organ.
- Utilize proper techniques to prevent air from entering the catheter system while removing and reinserting the pigtail wire. Improper techniques could result in an air embolism.
- The Instructions for Use supplied with all interventional devices to be used in conjunction with Baylis Medical Pigtail wires should be consulted for intended use, contraindications, warnings, precautions and potential complications related to such devices.
- When the pigtail wire is exposed to the vascular system, it should be manipulated while under high-resolution imaging guidance including fluoroscopy and/or echocardiography. Improper visualization of the guidewire may lead to misplacement, dissection, or perforation.
- Inspect the pigtail wire prior to use for coil separation, kinking, appropriate distal tip flexibility or breakage. If the pigtail wire is damaged or defective, do not use it. Using a damaged or defective pigtail wire may cause vasculature damage and/or compromise pigtail wire performance.
- Guidewire must always be advanced/withdrawn through a compatible catheter.
- The curve of the pigtail wire is not designed to be reshaped and should be handled with care to avoid damaging the distal tip or device.
- Laboratory staff and patients can undergo significant X-ray exposure during interventional procedures due to the

continuous usage of fluoroscopic imaging. The exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure.

- DO NOT attempt to deliver electrical signals or RF energy through device
- Failure to abide by the warnings in this labeling might result in damage to the device, which may necessitate intervention or result in serious adverse events

IV. PRECAUTIONS

- DO NOT allow the pigtail wire to remain in a prolapsed position as this may result in damage to the pigtail wire.
- DO NOT straighten the spiral distal tip of the wire more than 5 times during a procedure.
- Use care when removing the pigtail wire from the dispenser to reduce the possibility of pigtail wire damage.

V. ADVERSE EVENTS

Potential complications associated with the use of the pigtail wire include, but are not limited to:

- Vessel Perforation/Dissection/Trauma or Damage
- Vessel Spasm
- Hemorrhage
- Access Site Complications/Hematoma
- Thrombus/Thromboembolism
- Allergic reaction
- Vascular complication
- Cardiac tamponade
- Cardiac Perforation/Laceration
- Conduction disorder
- Embolism
- Additional Surgical Procedure
- Pericardial/pleural effusion
- Sepsis/Infection/Inflammation
- Foreign Body/Wire Fracture
- Hemolysis
- Hypovolemia
- Myocardial Ischemia and/or Infarction
- Stroke/Transient Ischemic Attack
- Vessel Occlusion
- Wire Entrapment/Entanglement
- Valve Complication

VI. INSPECTION PRIOR TO USE

Prior to use of the *ProTrack™* Pigtail Wire, the individual components should be carefully examined for damage or defects, as should all equipment used in the procedure. DO NOT use defective equipment. DO NOT reuse the device. The sterile packaging and kit containing devices should be visually inspected prior to use. DO NOT use the devices if the packaging has been damaged or compromised. Examine labelling to be intact and legible. Do not use if labelling is incomplete or illegible

VII. PREPARATIONS FOR USE

- Carefully read all instructions prior to use. Failure to do so may result in complications.
- Prepare other interventional devices according to the manufacturer's instructions.
- Remove the pigtail wire dispenser from the sterile package.
- Flush the dispenser with normal saline through the flushing port.
- Remove the pigtail wire insertion tool from the dispenser by removing the dispenser clip.
- Slowly advance the pigtail wire insertion tool over the pigtail wire's distal tip to straighten the spiral tip shape while holding the proximal end of the pigtail wire.
CAUTION: The tip of the pigtail wire is not designed to be reshaped and should be handled with care to avoid damaging the distal tip.
- Insert the pigtail wire insertion tool into the pigtail wire port of the catheter or sheath introducer. Carefully advance the pigtail wire's distal tip through the insertion tool into the pigtail wire lumen of the catheter. Remove the insertion tool by withdrawing it over the pigtail wire.

VIII. INSTRUCTIONS FOR USE

1. Carefully advance the pigtail wire through the pigtail lumen of the catheter until the distal end exits the catheter and recoils in the left atrium.
2. Confirm pigtail wire placement in two visual projections to ensure that the distal coil is fully extended from the tip of the catheter and in the correct position.
CAUTION: DO NOT advance or withdraw devices over the pigtail wire if the pigtail wire's spiral coil tip is not fully extended and visualized.
3. Utilize the pigtail wire to hold position while advancing catheters and sheaths over the pigtail wire into the desired location.
4. Withdraw the pigtail wire through the pigtail wire lumen of the catheter prior to attempting the final positioning of the catheter.

IX. SPECIAL STORAGE AND/OR HANDLING INSTRUCTIONS

Keep away from sunlight

X. CLEANING AND STERILIZATION INSTRUCTIONS

The pigtail wire is provided sterile and is intended for single-use only. DO NOT use if sterile barrier is damaged. If damage is found, device should not be used.

DO NOT clean or re-sterilize the pigtail wire. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

XI. DISPOSAL OF WASTE

Treat the used device(s) as biohazardous waste and dispose of in compliance with standard hospital procedures.

XII. CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION

If you have any problems with or questions about Baylis Medical Equipment contact our technical support personnel.



Baylis Medical Company Inc.

5959 Trans-Canada Highway
Montreal, Quebec, Canada, H4T 1A1
Phone: (514) 488-9801 or (800) 850-9801
Fax: (514) 488-7209
www.baylismedical.com

NOTES:

1. In order to return products you must have a return authorization number before shipping the products back to Baylis Medical Company. Product Return Instructions will be provided to you at this time.
2. Ensure that any product being returned to Baylis Medical has been cleaned, decontaminated and/or sterilized as indicated in the Product Return Instruction before returning it for warranted service. Baylis Medical will not accept any piece of used equipment that has not been properly cleaned or decontaminated as per the Product Return Instructions

XIII. LABELING AND SYMBOLS

	Catalogue Number		Contents
	Lot Number		Date of Manufacture
	EU Authorized Representative		Manufacturer
	Use By (Expiration Date)		Do Not Use if Packaging is Damaged
	See Instructions For Use		DO NOT re-sterilize
	Single Use Only (- DO NOT reuse)		Keep Away from Sunlight
	Sterilization with Ethylene Oxide Gas		Non-Pyrogenic
	Caution: Federal (U.S.A.) law restricts this device to sale by		EU Importer

	or on the order of a physician.		
	Quantity		Caution

XIV. LIMITED WARRANTIES – Disposables and Accessories

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date.

This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

XV. DISCLAIMER AND LIMITATION OF LIABILITY

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS

OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products:	The shelf life of the product
Accessory Products:	90 days from the shipment date

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