Instructions for Use

ProTrack™ Microcatheter

[English1]
 [Français2]
 [Deutsch4]
 [Italiano5]
 [Español6]
 [Português7]
 [Dutch9]

Baylis Medical Company Inc.
 5959 Trans-Canada Highway
 Montreal, Quebec, Canada, H4T 1A1
 Tel: (514) 488-9801 / (800) 850-9801
 Fax: (514) 488-7209
www.baylismedical.com

© Copyright [Baylis Medical Company Inc., 2007-2020]

EU Authorized Representative:

Quality First International
 20 Eversley Road, Bexhill-on-Sea, East Sussex
 TN40 1HE, United Kingdom
 Tel: +44-(20)-8-522-1937
 Fax: +44-(20)-8-522-1937

The Baylis Medical Company logo and Protrack are trademarks and/or registered trademarks of Baylis Medical Company Inc. in the United States of America and/or other countries.

English
 Carefully read all instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications. Baylis Medical Company relies on the physician to determine, assess and communicate to each individual patient all foreseeable risks of the procedure.

I. DEVICE DESCRIPTION

The ProTrack Microcatheter is designed to work as a standard infusion catheter. The catheter can also be used as an exchange guide for scaling-up thinner guidewires to thicker. All models have a removable luer-lock hub and are packaged with a removable stylet. These features permit the ProTrack Microcatheter to act as a guidewire. The critical dimensions of the catheter can be found on the device label.

The ProTrack Microcatheter is designed with a radiopaque distal tip marker that ensures precise tracking and placement of the catheter. The ProTrack Microcatheter also has a soft distal tip area that promotes easy tracking.

II. INDICATIONS FOR USE

The ProTrack Microcatheter is intended to be used as an exchange guide for converting guidewires or as a standard infusion wire.

III. CONTRAINDICATIONS

The ProTrack Microcatheter is not recommended for use with any condition that does not require an exchange guide for converting guidewires or a standard infusion wire.

IV. WARNINGS

- Only physicians with a thorough understanding of angiography and percutaneous interventional procedures should use this device.
- The ProTrack Microcatheter is supplied STERILE using an ethylene oxide process. Do not use if the package is damaged.
- The ProTrack Microcatheter is intended for single patient use only. Do not attempt to sterilize and reuse the catheter. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Failure to do so may result in patient complications.
- Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure.

V. PRECAUTIONS

- Do not attempt to use the ProTrack Microcatheter before thoroughly reading the accompanying Instructions for Use.
- Interventional procedures should be performed only by physicians thoroughly trained in the technique in a fully equipped catheterization laboratory.

- The sterile packaging should be visually inspected prior to use to detect any compromise. Ensure that the packaging has not been damaged. Do not use the equipment if the packaging has been compromised.
- Visually inspect the catheter to ensure there is no cracking or damage to the insulating material. Do not use the equipment if there is any damage.
- Do not use the ProTrack Microcatheter after the "Use By" date indicated on the label.
- Careful catheter manipulation must be performed to avoid cardiac damage, or tamponade. Catheter advancement should be done under fluoroscopic guidance. If resistance is encountered, DO NOT use excessive force to advance or withdraw the catheter.
- Do not bend or kink the catheter. Excessive bending or kinking of the catheter shaft may damage the integrity of the catheter and may cause patient injury. Care must be taken when handling the catheter.
- This device is not intended to be advanced through the vascular system without fluoroscopic guidance.
- Do not withdraw the ProTrack Microcatheter through a needle cannula.
- Do not exceed the maximum recommended pressure of 250 psi (1.72 MPa).
- Care must be taken to verify the dimensional compatibility of the ProTrack Microcatheter when used with other manufacturer's devices.
- Baylis Medical Company relies on the physician to determine, assess, and communicate to each individual patient all foreseeable risks of the ProTrack Microcatheter.
- Individual patient anatomy and physician technique may require procedural variations.

VI. ADVERSE EVENTS

Adverse events which may occur while using the ProTrack Microcatheter include:

Vessel perforation/dissection	Vascular thrombosis
Vessel spasm	Drug reactions
Hemorrhage	Allergic reaction to contrast medium
Hematoma	Arteriovenous fistula
Pain and Tenderness	Thromboembolic episodes
Death	Sepsis/Infection/Inflammation
Myocardial Infarction	Arrhythmias
Pericarditis	Pericardial/Pleural Effusion
Conduction/Heart Block	Cerebral Infarction
Air Embolus	Perforation and/or Cardiac Tamponade

VII. INSPECTION PRIOR TO USE

Prior to use, the ProTrack Microcatheter should be carefully examined for damage or defects, as should all equipment used in the procedure. Do not use defective equipment. Do not reuse the ProTrack Microcatheter.

VIII. EQUIPMENT REQUIRED

Intracardiac procedures should be performed in a specialized clinical setting equipped with a fluoroscopy unit, radiographic table, physiologic recorder, emergency equipment and instrumentation for gaining vascular access.

IX. SUGGESTED DIRECTIONS FOR USE

1. All instructions for equipment required should be carefully read, understood, and followed. Failure to do so may result in complications.
2. The ProTrack Microcatheter is supplied sterile. Remove the catheter with its attached luer-lock hub from the pouch using sterile technique.
3. Check luer-lock hub is firmly attached to catheter.
4. Thoroughly flush the catheter and auxiliary devices that may include guidewires and introducers with heparinized saline solution prior to use. Ensure luer-lock hub is tightened securely.
5. The ProTrack Microcatheter may be introduced into the vasculature through an introducer sheath, through a catheter, or over a wire. If resistance is encountered, DO NOT use excessive force to advance or withdraw the device. Determine the cause of resistance before proceeding. A compatible introducer device may assist introduction of the ProTrack Microcatheter through a hemostasis valve.
6. Use standard technique to position the catheter into the desired location under fluoroscopy. The radiopaque distal tip marker aids in fluoroscopic visualization during procedural placement.
7. Once the ProTrack Microcatheter has been negotiated into the desired position, infusion of fluid agents can begin in accordance with the manufacturer's instructions for dosage and duration. To inject/infuse fluid, withdraw the mandrill or guidewire and attach the infusion line or syringe to the luer-lock hub using standard technique.

NOTE: Before injecting fluid, always verify that the luer-lock hub is firmly attached to the ProTrack Microcatheter. DO NOT exceed the recommended maximum infusion pressure of 250 psi (1.72 Mpa).

To position the luer-lock hub, turn the proximal end counter-clockwise to loosen or clockwise to tighten. Proper placement of the luer-lock hub should be approximately 3 cm from the injectable proximal end of the ProTrack Microcatheter.

NOTE: The nominal flowrate of the catheter at 10PSI with saline solution is 10 ml/min and the nominal flowrate at 150PSI with saline solution is 16 ml/min.

X. CLEANING AND STERILIZATION INSTRUCTIONS

The ProTrack Microcatheter is intended for single use only and supplied sterile. Do not clean or re-sterilize the device. The device can be considered sterile only if the package is not opened or damaged prior to use.

XI. TROUBLESHOOTING

Bending and kinking of the catheter shaft may result in difficulty when maneuvering the catheter. Excessive bending and kinking may also cause patient injury. Care must be taken when handling the catheter. If bends, breaks, or kinks in the catheter are apparent, discard immediately.

XII. DISPOSAL OF WASTE

Treat used ProTrack Microcatheter as biohazard and dispose of in compliance with standard hospital procedures.

XIII. CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION

If you have any problems with or questions about Baylis Medical Equipment contact our technical support personnel.

Baylis Medical Company Inc.
5959 Trans-Canada Highway
Montreal, Quebec, Canada, H4T 1A1
Phone: (514) 488-9801 or (800) 850-9801
Fax: (514) 488-7209
www.baylismedical.com

NOTES:

- In order to return products you must have a return authorization number before shipping the products back to Baylis Medical Company. Product Return Instructions will be provided to you at this time.
- Ensure that any product being returned to Baylis Medical has been cleaned, decontaminated and/or sterilized as indicated in the Product Return Instruction before returning it for warranted service. Baylis Medical will not accept any piece of used equipment that has not been properly cleaned or decontaminated as per the Product Return Instructions.]

XIV. LABELING AND SYMBOLS

	Manufacturer
	EU Authorized Representative
	Sterile using ethylene oxide
	Use By
	Caution
	Consult Instructions for Use
	Model number
Rx ONLY	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
	Single Use – Do not reuse
	Do Not Resterilize
	Lot Number
	Do Not Use if Packaging is Damaged
	Keep Away From Sunlight
	Non-Pyrogenic
Max Guidewire O.D.	Maximum guidewire outside diameter that can be used with this device

XV. LIMITED WARRANTIES

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date. This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

DISCLAIMER AND EXCLUSION OF OTHER WARRANTIES

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date. This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

LIMITATION OF LIABILITY FOR DAMAGES

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products	The shelf life of the product
Accessory Products	90 days from the shipment date

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