Instructions for Use

RFP-100A Footswitch

Baylis Medical Company Inc.
5959 Trans-Canada Highway
Montreal, Quebec, Canada, H4T 1A1
Tel: (514) 488-9801/ (800) 850-9801 Fax: (514) 488-7209
www.baylismedical.com

EU Authorized Representative:
Quality First International
20 Eversley Road, Bexhill-on-Sea, East Sussex
TN40 1HE, United Kingdom
Tel: +44-(20)-8-522-1937 Fax: +44-(20)-8-522-1937

© Copyright Baylis Medical Company Inc., 2011-2019

The Baylis Medical logo is a trademark and/or a registered trademark of Baylis Medical Company Inc. in the United States of America and/or other countries.

English

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

I. DEVICE DESCRIPTION

The RFP-100A Footswitch (Footswitch) is used with the RFP-100A Baylis Medical Company Radiofrequency Puncture Generator (RFP-100A Generator) for RF energy delivery. The Footswitch is provided with redundant micro switches and is used as an alternative to the RF On/Off button on the RFP-100A Generator for RF energy delivery.

II. INDICATIONS FOR USE

The Footswitch is to be used with the RFP-100A Generator.

III. CONTRAINDICATIONS

The Footswitch is not recommended for uses other than the indicated use. The RFP-100A Generator Instructions for Use should be consulted for any other contraindications that may be associated with use of the RFP-100A Generator.

IV. WARNINGS

- Only physicians with a thorough understanding of radiofrequency powered puncture procedures should use this device.
- Do not alter this device in any way.
- The Footswitch must be used with the RFP-100A Generator only. Attempts to use it with other radiofrequency puncture generators can result in serious electrical injury of the patient and/or operator.

V. PRECAUTIONS

- Do not attempt to use the Footswitch before thoroughly reading the accompanying Instructions for Use and RFP-100A Generator Instructions for Use manual.
- Visually inspect the Footswitch prior to use. Do not use the Footswitch if there is any damage.
- Regularly inspect and maintain the Footswitch according to section VII. and X.

VI. ADVERSE EVENTS

The RFP-100A Generator Instructions for Use should be consulted for any adverse events.

VII. INSPECTION PRIOR TO USE

Prior to use all the individual components including the RFP-100A Generator and the Footswitch should be carefully examined for damage or defects, as should all equipment used in the procedure. Do not use defective equipment.

VIII. DIRECTIONS FOR USE

The Footswitch serves as an actuator for the connected device. When connected to the RFP-100A Generator, the Footswitch can be used as an alternative to the RF ON/OFF button on the RFP-100A Generator to activate or terminate RF energy delivery.

For Footswitch procedure and connections to the RFP-100A Generator please refer to the RFP-100A Generator Instructions for Use manual.

IX. CLEANING INSTRUCTIONS

The outer surface of the Footswitch may be cleaned with a mild soapy solution. DO NOT immerse the Footswitch in any liquid. Avoid caustic or abrasive cleaners. If disinfecting is required, 70% isopropyl alcohol or 5% solution of household bleach may be used to clean the outer surfaces. The Footswitch cannot be sterilized.

X. MAINTENANCE

Regular checks or maintenance is recommended according to ambient conditions and application frequency:

Check actuating elements for mechanical function. To perform this test, disconnect the footswitch from RFP-100A Generator. Press on the footswitch and release force. Confirm that the footswitch returns to its rest position (released state).

Check housing and connecting cable for damage and harmful dirt.

XI. CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION

If you have any problems with or questions about Baylis Medical Equipment contact our technical support personnel.

NOTES:

- In order to return products you must have a return authorization number before shipping the products back to Baylis Medical Company.
- Baylis Medical will not accept any piece of used equipment without a sterilization certificate.

XII. LABELING AND SYMBOLS

	Manufacturer		Consult Instructions for Use
	EU Authorized Representative		Model number
	Caution		Serial Number
	Only for EU member states: Use of this symbol indicates that the product must be disposed of in a way that complies with local and national regulations. For questions regarding recycling of this device please contact your distributor		
	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.		

XIII. Technical Data

Standards	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-2 IEC 60601-2-22 IEC 60529
Enclosure	Unbreakable, self extinguishing thermoplast.
Contact	Reed
Contact Material	Rhodium
Protection Class	IP X8 as per IEC 60529
Electrical Connection	3 m shielded cable (3xAWG 24C UL, GY) with LEMO-Connector 4 –pole.
Electrical Characteristics:	

Switching Voltage (UB)	Max. 25 V AC / 60 V DC
Low Current Switching (IB)	Max. 1A
Breaking Capacity (PB)	Max. 30VA
Mechanical lifetime	> 1Mio. Switching Cycles

XIV. Environmental Data

Environmental Conditions	Storage	Operating
Temperature	-20 °C to +70° C	-10 °C to +60 °C
Relative humidity	5 - 100%	10 - 100%
Air pressure	500hPa - 1060 hPa	700hPa - 1060 hPa

XV. LIMITED WARRANTY – Disposables and Accessories

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date. This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

DISCLAIMER AND LIMITATION OF LIABILITY

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products	The shelf life of the product
Accessory Products	90 days from the shipment date

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