

# Radiofrequency-Assisted Transseptal Access for Atrial Fibrillation Ablation Via a Superior Approach

Jackson J. Liang, DO, Aung Lin, MD, Sanghamitra Mohanty, MD, MS, Daniele Muser, MD, David F. Briceno, MD, J. David Burkhardt, MD, Gregory E. Supple, MD, David J. Callans, MD, Sanjay Dixit, MD, Rodney P. Horton, MD, Luigi Di Biase, MD, PHD, Francis E. Marchlinski, MD, Andrea Natale, MD, Pasquale Santangeli, MD, PHD

Liang et al. JACC: Clinical Electrophysiology, Jan 2020.

## INTRODUCTION

- ▶ This study reports outcomes of transseptal puncture and AF ablation from a superior approach in patients without access to the inferior vena cava (IVC) using dedicated radiofrequency (RF) tools.

## METHODS

- ▶ Retrospective analysis was performed on 15 patients undergoing RFA using the superior approach after an initial failed attempt using the femoral route.

### Transseptal puncture

- ▶ Ultrasound-guided superior access was obtained using the internal jugular, subclavian or axillary veins; radial arterial access was used for continuous blood pressure monitoring.
- ▶ SupraCross® Steerable Sheath (Baylis Medical) or Agilis™ EPI Steerable Introducer (Abbott) were advanced into the right atrium and deflected towards interatrial septum.
- ▶ A dedicated RF needle or wire (Baylis Medical) were used for transseptal puncture under intracardiac echocardiography (ICE) and fluoroscopy guidance.
- ▶ The steerable sheath was advanced into the left atrium over the pigtail SupraCross® RF Wire or ProTrack™ Pigtail Wire (Baylis Medical).

### Radiofrequency ablation

- ▶ Electroanatomic mapping (CARTO® 3 System, Biosense Webster) and wide antral PVI using a THERMACOOL® catheter (Biosense Webster) were performed.


## RESULTS

- ▶ Single (8 patients) or double (7 patients) transseptal access was obtained within  $16.1 \pm 4.8$  min.
- ▶ Mapping and ablation were performed successfully in 100% of cases within  $227.9 \pm 120.7$  min.

## DISCUSSION & CONCLUSIONS

- ▶ Superior transseptal access for AF ablation in patients with interrupted IVC can be achieved safely and effectively from a superior approach using dedicated RF transseptal devices.
  - May also provide a simpler strategy for delivery of endocardial LV leads

- ▶ Downward force on a standard transseptal needle and non-deflectable sheath from a superior approach can dislodge the contact site on the fossa ovalis.
- ▶ Use of a steerable sheath and stiff pigtail RF wire supported LA catheterization without difficulty.
- ▶ To better engage the fossa in two patients, the Agilis™ EPI sheath was exchanged for the SupraCross® sheath.
  - SupraCross® sheath had a smoother sheath-dilator transition than the Agilis™ sheath
  - SupraCross® sheath had a more flexible distal end that can be easily bent with the dilator in place for a tighter angle of deflection than the Agilis™ sheath (Figure 1)
- ▶ Technical considerations for performing AF ablation from the superior approach include:
  - ICE images are inverted compared to standard view; descending aorta can be used as a reference
  - Superior-anterior transseptal puncture improves catheter contact and stability in the left atrium
  - Left axillary vein access may improve operator ergonomics over right internal jugular vein



**Figure 1** Graphical reconstruction of sheath curvature presented in Liang et al using the (A) Agilis™ EPI sheath or (B) SupraCross® sheath with flexible dilator.