


**Baylis**  
MEDICAL

**TorFlex™**

# Transseptal Guiding Sheath and Dilator



# Torflex™

## TRANSSEPTAL GUIDING SHEATH AND DILATOR



Specifically designed to be used with the Toronto Transseptal Catheter for creating a transseptal puncture.

- Braided sheath construction provides superior torque control yet retains flexibility.
- Enhanced radiopaque tip maximizes visualization of sheath during manipulation.
- Tapered distal profile facilitates a smooth transition from dilator to sheath.
- Smooth inner lumen enhances ease of introduction and delivery of devices through sheath.

Catalog Number	Sheath French Size	Maximum Guidewire size Inches	Usable Length cm
TF8-38-62-S	8	0.038	62

© Copyright Baylis Medical Company Inc., 2007. Baylis Medical Company Inc. reserves the right to change specifications or to incorporate design changes without notice and without incurring any obligation relating to equipment previously manufactured or delivered. TorFlex™ and Baylis Medical Company Inc. logo are trademarks or registered trademarks of Baylis Medical Company Inc. in the United States of America and/or other countries.

CAUTION: Federal Law (USA) restricts the use of this device to or by the order of a physician.

CAR1006



Baylis Medical Company Inc.  
5959 Trans-Canada Highway  
Montreal, QC Canada H4T 1A1

Tel. : (514) 488-9801 / Fax : (514) 488-7209  
[www.baylismedical.com](http://www.baylismedical.com) / [info@baylismedical.com](mailto:info@baylismedical.com)