





**ProTrack™**  
Pigtail Wire



INSTANT ANCHOR. SECURE ACCESS.™

# ProTrack™ Pigtail Wire

Efficiently achieve and maintain left atrial access with a flexible spiral tip. Support large sheath exchange and cushion against the left atrial wall.




Radiopaque coil




### Save Time

Quicker introduction of devices into the LA.<sup>1</sup>  
Simply deploy the pigtail wire for instant, anchored access into the left atrium.



### Flexible Cushion

Cushion gently against the left atrial wall with a flexible, atraumatic tip.



### Exchange & Support

The stiff stainless steel body facilitates large sheath exchange and provides sheath stabilization.

## SPECIFICATIONS

Feature	Specifications
Body outer diameter	0.025 in
Overall length	175 cm, 230 cm
Functional exchange length	152 cm, 207 cm
Spring coil length	13 cm
Spiral coil diameter	4.1 cm
Shaft material	Stainless steel
Coil material	Stainless steel

## ORDERING INFORMATION

Overall Length	Product Code
175 cm	PTW-25-175
230 cm	PTW-25-230

<sup>1</sup>Buchner S, et al. J Interv Cardiol. doi: 10.1111/joic.12224

All trademarks are property of their respective owners. Patents Pending and/or issued. Caution: The law restricts this device to sale by or on the order of a physician. Indications, Contraindications, Warnings, and Instructions For Use can be found in the product labelling supplied with each device or at [www.baylismedical.com](http://www.baylismedical.com).

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries.  
This material not intended for use in France.

