




**ProTrack™**  
Microcatheter

FOR EXCELLENT TRACKABILITY.

# ProTrack™ Microcatheter


Use to exchange guidewires, infuse diagnostic or therapeutic agents, and facilitate access to distal tortuous vasculature over a guidewire.

## DESIGNED FOR VERSATILITY




### Excellent Trackability

Facilitate access to distal, tortuous vasculature over a guidewire.



### Easy Exchange

Exchange catheters effortlessly with a removeable hub.



### Radiopaque Tip

Flexible 10 cm distal end with a radiopaque tip marker visible under fluoroscopy.

## SPECIFICATIONS

Model	Outer Diameter Inches / French	Inner Diameter Inches	Length cm
CIC 35-145	0.035/2.7	0.022	145
CIC 38-145*	0.038/2.9	0.025	145

Maximum recommended pressure 250 psi

\*Compatible with Nykanen RF Wire

PRM-00092 EN J-1,2,3 V-2 © Copyright Baylis Medical Company Inc., 2013-2019. Baylis Medical Company Inc. reserves the right to change specifications or to incorporate design changes without notice and without incurring any obligation relating to equipment previously manufactured or delivered. ProTrack and the Baylis Medical logo are trademarks or registered trademarks of Baylis Medical Company Inc. in the United States of America and/or other countries. Before use, consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use. CAUTION: Federal Law (USA) restricts the use of this device to or by the order of a physician.

Products shown may not be approved in all jurisdictions.



# Baylis MEDICAL

Baylis Medical Company Inc.  
5959 Trans-Canada Highway  
Montreal, QC Canada H4T 1A1

Tel.: (514) 488-9801 / Fax: (514) 488-7209  
www.baylismedical.com / info@baylismedical.com