

The only
EXCHANGELESS
solution for access-to-delivery
of left heart ablation devices




**SINGLE SOLUTION
NO EXCHANGES**
deliver left heart ablation devices
with ease

VersaCross®
RF wire

VersaCross® dilator with
TRUform™ Shapeable
Technology

VersaCross® Steerable Sheath

3-in-1 RF wire facilitates left heart
access without exchanges:



VersaCross
RF Transseptal
Solution

Mechanical
needle

VersaCross
RF Transseptal Solution

1. Choose your **VersaCross®** RF Wire



Length: 180 cm, 230 cm
Diameter: 0.035"

2. Choose your **VersaCross®** sheath to
complete your solution



Available with Dilator Curves: D1, D0

VersaCross® Transseptal Sheath: Length: 63 cm, 81 cm Sheath Curves: 45°, 55°, 90° 8.5F	VersaCross® Steerable Sheath: Length: 71 cm Size: Sm, Md, Lg 8.5 F
---	--

PRM-00796 EN J-1,2,3 V-1 © Copyright Baylis Medical Company Inc., 2019-2021. Baylis Medical Company Inc. reserves the right to change specifications or to incorporate design changes without notice and without incurring any obligation relating to equipment previously manufactured or delivered. TRUform, OMNIVIZ, VersaCross, the VersaCross logo, Baylis Medical, and the Baylis Medical logo are trademarks and/or registered trademarks of Baylis Medical Company Inc. in the USA and/or other countries. Patents Pending and/or issued. CAUTION: Federal Law (USA) restricts the use of these devices to or by the order of a physician. Before use, consult product labels and Instructions for Use for Indications for Use, Contraindications, Warnings, Precautions, Adverse Events and Directions for Use.

Products shown may not be approved in all jurisdictions.



Baylis
MEDICAL

Baylis Medical Company Inc.
5959 Trans-Canada Highway
Montreal, QC Canada H4T 1A1

Tel.: (514) 488-9801 / Fax.: (514) 488-7209
www.baylismedical.com / info@baylismedical.com

Baylis
MEDICAL

VersaCross
RF Transseptal Solution



**EP WORKFLOW
ADVANTAGES**



Mechanical needle