


# FAST TRACK

to your therapy delivery  
in a SINGLE SOLUTION



## VersaCross Connect™ LAAC ACCESS SOLUTION

DESIGNED for use with  
**WATCHMAN™**  
Access Sheaths<sup>‡</sup>


**3-IN-1  
RF WIRE**  
for ZERO EXCHANGE delivery  
of WATCHMAN™ Access Sheath<sup>‡</sup>



Smooth transition between  
VersaCross Connect™ Dilator and  
WATCHMAN™ Access Sheath<sup>‡</sup>



Snap-lock designed for  
WATCHMAN™ Access Sheaths<sup>‡</sup>



1. Choose your VersaCross® RF Wire  
Pigtail  
Length: 180 cm, 230 cm  
Diameter: 0.035"
2. Choose your VersaCross Connect™ Transseptal Dilator  
Dilator with TRUform™ technology  
D0  
D1  
Inner Diameter: 0.035"  
Outer Diameter: 12F  
Useable Length: 85cm

### Solution with Pigtail Wire

Part No.	Dilator curve	Wire length (cm)
VXAK0001	D0	180
VXAK0003	D1	180
VXAK0005	D0	230
VXAK0007	D1	230

<sup>‡</sup> The VersaCross Connect™ Transseptal Dilator is for use with a 12F ID WATCHMAN™ Access Sheath that is 75 cm in length, specifically: WATCHMAN™ Access System [Models: M635TU40060, M635TU10060, M635TU20060]; WATCHMAN™ TruSeal™ Access System [Models: M635TU70010, M635TU70040, M635TU70020]; WATCHMAN FXD Curve™ Access System [Models: M635TU80010, M635TU80020].

VersaCross Connect™ LAAC Access Solution includes VersaCross Connect™ Transseptal Dilator and VersaCross® RF Wire.  
PRM-00907 EN J-1,2 V-2 © Copyright Baylis Medical Company Inc., 2022. Baylis Medical Company Inc. reserves the right to change specifications or to incorporate design changes without notice and without incurring any obligation relating to equipment previously manufactured or delivered. TRUform, OMNIviz, VersaCross, VersaCross Connect, Baylis Medical and the Baylis Medical logo are trademarks and/or registered trademarks of Baylis Medical Company Inc. in the USA and/or other countries. Other trademarks property of their respective owners. Patents Pending and/or issued. Before use, consult product labels and Instructions for Use for Indications for Use, Contraindications, Warnings, Precautions, Adverse Events and Directions for Use.

Products shown may not be approved in all jurisdictions.



Baylis Medical Company Inc.  
5959 Trans-Canada Highway  
Montreal, QC Canada H4T 1A1  
Tel.: (514) 488-9801 | Fax.: (514) 488-7209  
www.baylismedical.com | info@baylismedical.com



## ZERO EXCHANGE WATCHMAN™ SHEATH DELIVERY

Connects VersaCross®\* technology  
with WATCHMAN™<sup>‡</sup> sheaths for  
ZERO EXCHANGE left heart access

\* VersaCross Connect™ LAAC Access Solution  
<sup>‡</sup> WATCHMAN™ Access System



\*\* WATCHMAN™ Access Sheaths and devices sold separately

# ZERO EXCHANGE

WATCHMAN™ Access Sheath Delivery†

## EXCHANGELESS 3-in-1 RF Wire



1 Start at the SVC



2 RF transseptal puncture




3 Deliver WATCHMAN™ Sheath

INSTANTLY DEPLOY  
**0.035" WIRE**

to deliver  
**WATCHMAN™ Access Sheath†**  
with confidence

## DELIVER THERAPY ON TARGET



**TRUform™**  
Shapeable Technology



**Precision RF**  
Puncture Technology



**OMNIviz™**  
Technology

Reliably locate your **VersaCross Connect™** solution on fluoroscopy, ultrasound and visibly with positional markers

the right **CURVE** at the right **TIME**



**TRANSSEPTAL CURVE**

**VersaCross Connect™**  
Transseptal Dilator with **TRUform™** shapeable technology to optimize transseptal approach



**THERAPY DELIVERY CURVE**

**WATCHMAN™ Access Sheath** retains its curve for therapy delivery

 **VersaCross**  
Connect Access Solution